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Collaborator: Nippon Becton Dickinson Company, Ltd.

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Flow cytometer is a device that can semi-quantitatively measure the expression of individual cell molecules. There are two types of devices, an analyzer and a cell sorter. Antibodies recognized by fluorescent dyes are often used to identify the expression of molecules. It has also been attempted to measure fluorescent substances other than labeled antibodies.

As a practical training, we detected with FACS using beads for flow cytometer setup.

* FACS (fluorescence activated cell sorter) is a trademark of Becton, Dickinson and company.

Practice content:

Beads for flow cytometer setup are detected with FACS.

We will analyze using three instrument, analyzer BD LSRFortessaX-20 and BD FACSCanto II, and cell sorter BD FACSAria Fusion.

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Explanation of instrument

Data acquisition and analysis

Practice will be done in 3 groups and 2 instruments in turn.

Group A; LSRFortessa (17:00–17:40) → AriaFusion (17:40–18:20)

Group B; Canto II (17:00–17:40) → Aria (17:40–18:20)

Group C; Aria (17:00–17:40) → Canto II (17:40–18:20)

The aim for today's practice is to learn about the settings and analysis for FACS.

We want you to know that FACS is a useful and easy instrument.

■ Practice

Samples to measure;

Samples are prepared in the tubes as shown below.

Adjust the voltage/compensation for settings.

Analyze the percentage of PE positive Beads.

Tube NO.	1	2	3	4	5
label	Unlabeled	Unlabeled FITC	Unlabeled PE	Unlabeled APC	Unlabeled FITC PE APC
Unlabeled	○	○	○	○	○
FITC		○			○
PE			○		○
APC				○	○