

How to Use Flow Cytometer and Cell Sorter (Lecture in English)

幾尾 真理子、大家 吉泰 (ウォーターズバイオサイエンス)
Mariko Ikuo, Yoshiyasu Ooie (Waters Biosciences)

フローサイトメーターは個々の細胞の分子の発現を半定量的に測定できる装置である。装置は解析専用のアナライザーと解析分取が可能なセルソーターの 2 機種があり、分子の発現の同定には蛍光色素で認識された抗体がよく用いられる。更に標識抗体以外の蛍光物質を測定する事も試みられている。今回は実習として、機器更正用ビーズを用いて、FACS にて検出する。

FACS の基本原理を紹介し、FACS を用いた実験手技については具体例を挙げながら説明します。

Flow cytometer is a device that can semi-quantitatively measure the expression of individual cell molecules. There are two types of devices, an analyzer and a cell sorter. Antibodies recognized by fluorescent dyes are often used to identify the expression of molecules. It has also been attempted to measure fluorescent substances other than labeled antibodies.

As a practical training, we detected with FACS using beads for flow cytometer setup.

* FACS (fluorescence activated cell sorter) is a trademark of Becton, Dickinson and company. I will introduce the basic principle of FACS and explain the experimental technique using FACS with concrete examples.

Waters Biosciences (formerly BD Biosciences)



For Research Use Only. Not for use in diagnostic or therapeutic procedures.
APCC® (patent: 5,714,386)
Aqua Pure® Catalyst and Texas Tech Blue® are registered trademarks and Pacific Blue™ and Pacific Orange™ are trademarks of Life Technologies Corporation.
Cyt is a trademark of Becton, Inc.
Cyt is a trademark of GE Healthcare. Cyt dyes are subject to proprietary rights of GE Healthcare and Carnegie Mellon University and are made and sold under license from GE Healthcare only for research and in-vitro diagnostic use. Any other use requires a commercial sublicense from GE Healthcare, 800 Centre Avenue, Pittsburgh, PA 15265-1327, USA.
Cyto is a registered trademark of eBioscience, Inc.
© 2016 Life Technologies Corporation. All rights reserved. All trademarks are property of Becton, Dickinson and Company. © 2016 BD
23-167-2901

NO RESTRICTION

BD Biosciences is now a part of Waters Corporation

This business has been acquired by Waters Corporation ("Waters"), Becton, Dickinson and Company or one of its affiliates or subsidiaries ("BD") remains the legal manufacturer of Biosciences and Diagnostic Solutions products until all required regulatory transfers are completed. During this interim period, BD maintains full responsibility for all regulatory obligations of the legal manufacturer. Product information is supplied under BD's regulatory authority. To learn more about the relationship between Waters and BD, please visit <https://www.bd.com/our-relationship>. BD is a trademark of Becton, Dickinson and Company. Waters is a trademark of Waters Corporation or its affiliates. All other marks are the property of their respective owners. © 2026 Waters Corporation.

BD is a trademark of Becton, Dickinson and Company. Waters is a trademark of Waters Corporation or its affiliates. All other marks are the property of their respective owners. © 2026 Waters Corporation.

BD is a trademark of Becton, Dickinson and Company. Waters is a trademark of Waters Corporation or its affiliates. All other marks are the property of their respective owners. © 2026 Waters Corporation.




Introduction to Flow Cytometry (FCM)

1. Basic of FCM
2. Application for FCM



3

Technique to analyze the character of cells or particle by measuring the scatter light and fluorescence exposed to a laser running in a stream



The diagram illustrates the workflow of a Cell Sorter, which is part of a larger system including a Cell Analyzer. The process is divided into five numbered steps:

- Sample flow**: A vertical tube where the sample is being processed.
- Laser hits the cell**: A laser beam (represented by a blue arrow) passes through the sample flow.
- Detect Scatter light and fluorescence**: Multiple colored arrows (red, yellow, green, blue, purple) represent the detection of different light signals from the cells.
- Display data**: A screenshot of a flow cytometry analysis software showing two plots: a scatter plot of $\log_{10}(\text{SSC-A})$ vs $\log_{10}(\text{FSC-A})$ and a histogram of $\log_{10}(\text{FSC-A})$. Below the plots is a table of statistics for various populations.
- Collect specific cells**: An orange arrow points from the selected population in the histogram to a collection tube containing sorted cells.

Cell Analyzer

Display data

Statistics table from the Cell Analyzer display:

Population	%Area	%Count	%Area	%Count
1-100	29.0	29.0	29.0	29.0
2-100	4.88	4.87	4.88	4.88
3-100	5.474	5.47	5.47	5.47
4-100	7.95	7.95	7.95	7.95
5-100	12.0	12.0	12.0	12.0
6-100	12.0	12.0	12.0	12.0



FOM15: • • •

- **In the late 1960s**
Leonard Herzenberg, PhD, professor at Stanford University developed Fluorescence Activated Cell Sorter (FACS)
- **1973**
Commercialized by US Becton-Dickinson (FACS I)
- **New products are released every few years**
Improving laser and data processing system, making it possible to analyze more items much faster



Instrument	Maximum number of laser equipment	Number of fluorescent color	Processing speed (Cells/sec)
Early FACS	1	2	-
Present FACS	6	50 or more	70000




Central Research Laboratory : Shiga University of Medical Science

Cell Sorter



Cell Analyzer




More Colors
&
Higher Performance




7
NO RESTRICTED

Character of FCM


FCM



- ✓ High Throughput
- ✓ Multiparameter
- ✓ Quantitative data
- ✓ Cell sorting
- ✓ Objectivity, Reproducibility
- ✓ Imaging (Imaging Cytometer)



Fluorescence microscope




- ✓ Image (3D)
- ✓ Intracellular distribution
- ✓ Time-lapse



8
NO RESTRICTED


FCM analysis targets

- Single cells in suspension
- Cell size 0.5 - 40 μm



9
NO RESTRICTED

Detection of FCM



Laminar Flow formation / Laser / Detection

Formation of narrow and laminar flow, enables cells to run one by one and passes through the laser


- Detects scatter and fluorescence
- Process few thousand cells/sec



Laminar Flow formation / Laser / Detection


10
NO RESTRICTED

Data Display : Dot Plot and Histogram



11
NO RESTRICTED

Character of scatter light




- Forward Scatter (FSC)
Indicates the size of cells, surface area
- Side Scatter (SSC)
Indicates the complexity of internal structure, cellular granule




12
NO RESTRICTED


Scatter light signal




Scatter light analysis of blood cells




Character of fluorescence




Fluorescent staining with monoclonal antibody




Fluorescent signal




Fluorescent signal analysis



Fluorescence detector of FCM




Compensation: Separating fluorochromes with spillover



Sorting Mechanism

1. Detect cells by scatter and fluorescence signals
2. Identify target cells using gating criteria
3. Generate droplets containing individual cells
4. Charge target droplets at droplet breakoff
5. Deflect charged droplets into collection tubes using an electric field
6. Collect purified cell populations for downstream analysis




Today's Contents

Introduction to Flow Cytometry (FCM)

1. Basic of FCM
2. Application for FCM

Applications for FCM

- ① Cell surface antigen
- ② Intracellular antigen
- ③ Fluorescent protein
- ④ Cell cycle
- ⑤ Cell proliferation
- ⑥ Apoptosis
- ⑦ Cytokine




① Analysis of cell surface antigen

Major cell surface antigen


- ✓ Cellular membrane composing protein
- ✓ receptors
- ✓ Adhesion molecules
- ✓ Transporters

Cluster of Differentiation (CD-classify)

International classification of monoclonal antibodies used for cell surface antigen



Analysis of cell surface antigen



25

NO RESTRICTED

② Analysis of intracellular antigen


- Intracellular antigens
 - Cytokine
 - Transcription factor
 - Cytoskeleton

Cells fixed and permeabilized makes antibody and reagents possible to access intracellular

- Fixation reagent
 - Ethanol
 - Methanol
 - Paraformaldehyde

- Permeabilization reagent
 - Saponin
 - Triton-X

- Nucleic acid staining reagent
 - PI
 - 7AAD




25

NO RESTRICTED

Cytokine producing cells analysis

- Staining with cell surface antigen, enables phenotype analysis of cytokine producing cells
- Possible to detect cytokine production even from rare cells




27

NO RESTRICTED

Protein phosphorylation analysis in single cell


- Protein phosphorylation in single cell level
- Combining the staining of cell membrane markers, allows protein phosphorylation analysis of specific cells in heterogeneous cell population
- Easy, short time and less cell sample compared to Western Blot analysis



28

NO RESTRICTED

Protein phosphorylation analysis in single cell




29

NO RESTRICTED

③ Expression of fluorescent protein


- Purpose
 - Transgene expression
 - Transplanted cells survival
 - Interaction between transgene
- Fluorescent proteins
 - GFP series fluorescent protein
 - GFP, YFP, CFP, RFP, BFP etc.
 - Fruits series fluorescent protein
 - mCherry, mPlum, mStrawberry, mBanana etc.
 - CoralHue series fluorescent protein
 - Midoriishi-Cyan, Kusabira-Orange, Azami-Green, Keima-Red, Kaede, Doronpa-Green etc.



30

NO RESTRICTED


GFP / mCherry Expressing cell analysis



31

NO RESTRICTED


GFP / YFP / CFP Expressing cell analysis



32

NO RESTRICTED

④ Cell cycle analysis



Cell cycle	Marker/reagent	mechanism
All	PI, 7-AAD	interacts with double stranded nucleic acid
S phase	BrdU/anti-BrdU Ab	BrdU replaces Thymidine(T) on DNA synthesis, BrdU is detected with Ab
G0/G1 phase	Anti-CyclinAb, anti-Rb Ab markers for cell cycle	Amount changes in the phases
M phase	Anti Phospho-HistoneH3 antibody	changes in phosphorylation




33

NO RESTRICTED

④ Cell cycle analysis : PI analysis

- PI staining cell cycle analysis
CycleTEST™ PLUS DNA Reagent Kit (Cat#340242)



Hela cells cultured 48hr in DMEM with 0.5% FBS were treated with reagents.
After 16hr culture, cell cycle analysis was performed




34

NO RESTRICTED

④ Cell cycle analysis : BrdU / 7AAD analysis

BD Pharmingen™ FITC BrdU Flow Kit (Cat#559619)




35

NO RESTRICTED

⑤ Cell proliferation assay



BD Biosciences offers BD Horizon™ Violet Proliferation Dye 450 (VPD450) and BD Horizon™ CFSE for the detection of cell proliferation. Once the dye is inside the cell, esterases cleave off the ester group to convert the dye into a fluorescent product and trap it inside the cell. With each replication event the amount of dye in the cell is decreased, leading to a characteristic pattern.



36

NO RESTRICTED

⑤ Cell proliferation assay: VPD450 assay



37

NO RESTRICTED

⑥ Apoptosis analysis

Feature Measured	Assays	Key Features
Plasma Membrane Alterations (Phosphatidylserine Exposure)	- Annexin V binding assay - Single conjugates - Annexin V Kits	- Detects early apoptosis markers - Quick and easy - Flow cytometry or immunofluorescence application
Mitochondrial Changes	BD MitoScreen Kit	Fast, easy, single cell resolution by flow cytometry or fluorescent microscopy
Caspase Activation	- Caspase Activity Assay Kits and Reagents - Active Caspase-3 immunoassays (ELISA)	- Quick and easy, uses spectrofluorometry - ELISA application
DNA Fragmentation	- APO-BrdU TUNEL Assay - APO-DIRECT TUNEL Assay	Works adherent cells, single cell resolution in conjunction with cell cycle analysis by flow cytometry




Many ways to analyze apoptosis

38

NO RESTRICTED

Apoptosis analysis: Annexin V assay for Surface Phosphatidylserine detection




Phosphatidylserine (PS), which is normally located on the cytoplasmic face of the plasma membrane, translocates to the outer leaflet of the plasma membrane during apoptosis and can be detected by flow cytometry and cell imaging through binding to fluorochrome-labeled annexin V when calcium is present.



39

NO RESTRICTED


Apoptosis analysis: Annexin V assay



40

NO RESTRICTED

Apoptosis analysis: Active Cleaved Caspases and Cleaved poly-ADP ribose polymerase (PARP) detection



Caspases are activated upon during the earliest stages of apoptosis. Active caspases can then cleave many proteins, including PARP. This leads to the loss of cellular structure and function, and ultimately results in cell death.




Detects activated Caspase3 with specific antibody

41

NO RESTRICTED

Apoptosis analysis: Apo-Direct to detect DNA damage



One of the later steps in apoptosis is DNA fragmentation. A method which is often used to detect fragmented DNA utilizes a reaction catalyzed by exogenous TdT, often referred to as "end-labeling" or "TUNEL" (terminal deoxynucleotidyltransferase dUTP nick end labeling).

The APO-DIRECT™ assay is a single-step method for labeling DNA breaks with FITC or APC-dUTP.



Detects DNA damage in apoptosis cells


42

NO RESTRICTED


⑦ Cytokine measurement

method	object	measurement	instrument
ELISA (OptEIA)	Cytokine concentration in solution	Antibody and enzyme reaction in 96 well plate to color	Plate reader
Intracellular Cytokine	Cells producing cytokine	Stain cytokine with Ab within cells	FCM
ELISPOT	Cells secreting cytokine	Antibody and enzyme reaction in 96 well plate to color	Microscope
Cytometric Beads Array (CBA)	Cytokine concentration in solution	Use special beads to measure cytokine	FCM

Cytometric Beads Array (BD™ CBA)



Cytometric Beads Array (BD™ CBA)



	BD™ CBA	ELISA
Multiplexing Capability	Up to 30 proteins /sample	one protein /sample
Simple procedure and Time Efficiency	Wash once	Wash more than 10 times
Sample amount	50 µl or less	About 100 µl
Format	From single tube	96 well plate
Required Instrument	Flow cytometer	Plate reader

Cytometric Beads Array (BD™ CBA)



Capture beads for each different target protein are distinguished by unique fluorescence intensity of beads. The PE intensity of each capture bead population reflects the target protein concentration. Standard curves are plotted using titrated standard protein samples.

Searching Waters Biosciences products



Enter Keyword Here!

*Access our web site →
<https://www.bdbiosciences.com/ja/jp>
<https://www.bdbiosciences.com/en/us>

Jumps to list of related products



- ✓ Click catalog number for details
- ✓ Available to download datasheet
- ✓ Prices are shown, useful for ordering

Guide for multicolor analysis, posters, protocol, etc., are available to download

Contact Details

Flow cytometry operator trainings, application trainings and FlowJo™ training

Training takes place at our lab in Tokyo or at your lab.
Please contact for details.

E-mail : customer_training@bd.com

For question about instruments and reagents

Please contact application support
for any questions about our instruments and reagents

E-mail : tech_cell@bd.com

For machinery problems and maintenance

Instrument support

<https://bdbqcodeservice.powerappsportals.com/>