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


# Spatial Transcriptomics is Hot, Rapidly emerged field

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
Visium(55um)

VS

In his view, researchers will want single-cell resolution, which scRNA-seq provides but tissue imaging, as of yet, does not. Resolution matters because in a region of interest tissues can have a mix of cell types but "a completely wrong picture" can emerge

# PhenoCycler is the leading Spatial Phenotyping solution for Global Cell Atlas Efforts

Global program to create cell atlas


Top tier institute participation

Pacific Northwest  
NATIONAL LABORATORY

UC San Diego

Nearly 100 datasets contributed to public domain  
Over 50 publications

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

**Autofluorescence  
Reduction:**  
**LED + Incubation**  
**4.5% H<sub>2</sub>O<sub>2</sub> +**  
**20mM NaOH in PBS**  
**2x 45 min**


# Quantitative Evaluation of Autofluorescence Reduction

FFPE Human Brain: Two serial sections

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

Human FFPE			Human FF	Mouse FF
Anti-Hu ATM	Anti-Hu CD107a	Anti-Hu LAG3	Anti-Hu CD2	Anti-Mu CD3
Anti-Hu b-Catenin1	Anti-Hu CD141	Anti-Hu LIF	Anti-Hu CD3	Anti-Mu CD4
Anti-Hu Bcl-2	Anti-Hu CD163	Anti-Hu Mac2/Galectin-3	Anti-Hu CD4	Anti-Mu CD5
Anti-Hu Beta-actin	Anti-Hu Collagen IV	Anti-Hu MPO	Anti-Hu CD8	Anti-Mu CD8a
'Anti-Hu Beta-actin	Anti-Hu E-cadherin	Anti-Hu Pan-Cytokeratin	Anti-Hu CD11c	Anti-Mu CD11b
Anti-Hu Caveolin	Anti-Hu EpCAM	'Anti-Hu Pan-Cytokeratin	Anti-Hu CD19	Anti-Mu CD11c
Anti-Hu CD3e	Anti-Hu ER	Anti-Hu PCNA	Anti-Hu CD21	Anti-Mu CD19
Anti-Hu CD4	Anti-Hu FOXP3	Anti-Hu PD-1	Anti-Hu CD31	Anti-Mu CD21/35
Anti-Hu CD8	Anti-Hu GATA3	Anti-Hu PD-L1	Anti-Hu CD34	Anti-Mu CD24
Anti-Hu CD11c	Anti-Hu GP100	Anti-Hu Podoplanin	Anti-Hu CD38	Anti-Mu CD31
Anti-Hu CD14	Anti-Hu Granzyme B	Anti-Hu SMA	Anti-Hu CD45	Anti-Mu CD38
Anti-Hu CD20	Anti-Hu HIF1A	Anti-Hu SOX2	Anti-Hu CD45RO	Anti-Hu/Mu CD44
Anti-Hu CD21	Anti-Hu Histone H3 Phospho	Anti-Hu T-bet/TBX21	Anti-Hu/Mu CD49f	Anti-Mu CD45
Anti-Hu CD31	Anti-Hu HLA-A	Anti-Hu TCF-1	Anti-Hu CD69	Anti-Mu CD45R/B220
Anti-Hu CD34	Anti-Hu HLA-DR	Anti-Hu TFAM	Anti-Hu CD90	Anti-Hu/Mu CD49f
Anti-Hu CD38	'Anti-Hu HLA-DR	Anti-Hu TOX	Anti-Hu CD104	Anti-Mu CD71
Anti-Hu CD39	Anti-Hu HLA-E	Anti-Hu TP63	Anti-Hu CD138	Anti-Mu CD90.2
Anti-Hu CD40	Anti-Hu ICOS	Anti-Hu Vimentin	Anti-Hu CD278	Anti-Mu CD169
Anti-Hu CD44	Anti-Hu IDO1	Anti-Hu VISTA	Anti-Hu HLA-DR	Anti-Mu IgD
Anti-Hu CD45	Anti-Hu IFNG		Anti-Hu/Mu Ki67	Anti-Mu IgM
Anti-Hu CD45RO	Anti-Hu iNOS		Anti-Hu Pan-Cytokeratin	Anti-Hu/Mu Ki67
Anti-Hu CD56	Anti-Hu Keratin 5		Anti-Hu PD-1	Anti-Mu Ly6g
Anti-Hu CD66	Anti-Hu Keratin 8/18		Anti-Hu Podoplanin	Anti-Mu MHC II
Anti-Hu CD68	Anti-Hu Keratin 14			Anti-Mu TCRb
Anti-Hu CD79a	Anti-Hu/Mu Ki67			Anti-Mu Ter119

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
More flexible in cell  
segmentation options →

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## Public open source

bioRxiv

### Evaluation of cell segmentation methods without reference segmentations

Haoran Chen, Robert F. Murphy

doi: <https://doi.org/10.1101/2021.09.17.460800>

Table 3 Segmentation methods evaluated

Method	Inputs					Output
	Cytoplasm	Cell Membrane	Nucleus	Dimension	Requires scale	Nuclear mask
DeepCell 0.6.0		X	X	2D	No	Yes
DeepCell 0.6.0	X		X	2D	No	Yes
DeepCell 0.9.0 (new)		X	X	2D	No	Yes
DeepCell 0.9.0 (new)	X		X	2D	No	Yes
Cellpose 0.0.3.1	X		X	2D/3D	No	Yes
Cellpose 0.6.1 (new)	X		X	2D/3D	No	Yes
CellSegm		X		2D/3D	Yes	No
CellIX		X		2D	Yes	No
CellProfiler		X	X	2D	Yes	No
AICS (classic)		X	X	2D/3D	Yes	No
Voronoi			X	2D	No	Yes

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


Phenotype Assignment	Population	Positive Markers (All CD45+)
Stromal Cells	3414	CD90+, CD9+, HLA-DR+, Cytokeratin+
Stromal Cells in connective regions	807	CD31+, CD90+, Collagen IV+
T <sub>H</sub> Cells	328	CD3+, CD4+, CD7+, CD278+, CD279+
Activated T <sub>H</sub> Cells	1459	CD3+, CD4+, CD40+, CD278+, CD279+
Activated Helper T cells	11866	CD3+, CD4+, CD7(+)
T Cells near connective regions	222	CD3+, CD29+, CD90(+), Collagen IV(+)
Proliferating Th-1 cells	2411	CD4+, Ki67+
Dendritic Cells	135	CD4+, HLA-DR+
Dendritic Cells	366	CD4+, CD11c+
Dendritic Cells	429	CD11c+, HLA-DR+
NK cells	126	CD3-, CD7+
Plasma Cells	2180	CD31+, CD38+, CD40+
Lymphatic Vasculature Cells	1336	CD9+, Podoplanin+
Mature, proliferating B cells near follicular matrix	960	CD19+, CD22+, CD21+, CD90+, Ki67+
Mature B cells	32840	CD19+, CD22+

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Cell Proximity Heat-map


More interacting  
Less interacting

- Plasma cells
- Plasma cells near epithelial cells
- DCs near B cells
- NKs
- T cells interacting with DCs
- Vascular endothelial cells (vWF/CD146)
- Combined T cell population
- Epithelial cells (CD63/Pan-cytokeratin)
- Vascular endothelial cells (CD146/CD36)
- Vascular endothelial cells (CD146/vWF/CD36)
- B cells
- DCs interacting with T cells
- Epithelial cells (CD63/vimentin/Pan-cytokeratin)
- Epithelial cells (CD63/CD4/Pan-cytokeratin)
- T cells interacting with DCs
- Helper T cells
- Combined T cell population
- T cells interacting with DCs
- TH1 (CD279+)
- Lymphatic vasculature
- NKs
- Plasma cells
- Plasma cells near epithelial cells
- Vascular endothelial cells (vWF/CD146)
- Vascular endothelial cells (CD146/CD36)
- Filtered

# Cellular Neighborhood解析 (CytoMAP)



<https://www.akoyabio.com/webinar/codex-with-seurat-a-high-dimensional-approach-to-unlocking-spatial-insights/>

## Cellular Neighborhood解析


Cluster spatial windows to find common patterns of proximity, i.e., cell neighborhoods

- ● ● ● ● ● ● ● ● ● Cell 1 neighbors (window 1)
- ● ● ● ● ● ● ● ● ● Cell 2 neighbors (window 2)
- 
- 
- 

## Cellular Neighborhood Interaction

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




Cell

### Coordinated Cellular Neighborhoods Orchestrate Antitumoral Immunity at the Colorectal Cancer Invasive Front

Cell. 2020 Oct 29;183(3):838.



Resource

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CODEX multiplexed tissue imaging with DNA-conjugated antibodies

Sarah Black<sup>1,2,6</sup>, Darci Phillips<sup>1,2,6</sup>, John W. Hickey<sup>1,2,6</sup>, Julia Kennedy-Darling<sup>1,4</sup>, Vishal G. Venkataraman<sup>1,2</sup>, Nikolay Samusik<sup>1,5</sup>, Yuri Goltsev<sup>1,2</sup>, Christian M. Schürch<sup>1,2,3,7,8</sup> and Garry P. Nolan<sup>1,2,7,8</sup>

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

Model of differences in the iTME between CLR and DII patients with respect to CN organization, cellular function, and inter-CN communication

# New Biomarkers are Needed to Better Predict Patient Response

ARTICLE

<https://doi.org/10.1038/s41467-021-26974-6>


OPEN

Check for updates

## Immune cell topography predicts response to PD-1 blockade in cutaneous T cell lymphoma

Darci Phillips<sup>1,2,3,9</sup>, Magdalena Matusiak<sup>3,9</sup>, Belén Rivero Gutierrez<sup>3</sup>, Salil S. Bhate<sup>1,3,4</sup>, Graham L. Barlow<sup>1,3</sup>, Sizun Jiang<sup>1,3,5</sup>, Janos Demeter<sup>1</sup>, Kimberly S. Smythe<sup>6</sup>, Robert H. Pierce<sup>6</sup>, Steven P. Fling<sup>6</sup>, Nirasha Ramchurren<sup>6</sup>, Martin A. Cheever<sup>6</sup>, Yury Goltsev<sup>1,3</sup>, Robert B. West<sup>3</sup>, Michael S. Khodadoust<sup>7,10</sup>, Youn H. Kim<sup>2,7,10</sup>, Christian M. Schürch<sup>1,3,8,10</sup> & Garry P. Nolan<sup>1,3,10</sup>


## Conventional Biomarkers Fail to Predict Outcomes in Lymphoma

# Context Matters: Spatial Neighborhoods Predict Response



## Discovery Uncovers Predictive Potential of Spatial Neighborhoods

### CITN-10 Phase II Clinical Trial (Cutaneous T-Cell Lymphoma)

- 56-plex IO Panel
- Multiple spatial neighborhoods differentiate responders vs. non-responders!

Spatial Neighborhood Map Responder

# Spatial Neighborhoods in Today's Literature

## Rapid Adoption of SNs as Potential Biomarkers Across the Biomedical Literature

**Cell** Resource

### Coordinated Cellular Neighborhoods Orchestrate Antitumoral Immunity at the Colorectal Cancer Invasive Front

Graphical Abstract

**Technology**  
CODEX + Optimization for Formalin Fixed Paraffin Embedded Tissue (FFPE)

**Data Collection**  
Colorectal Cancer (High Risk, Low Risk) → FFPE → Tissue Microarray

**Computation**  
Cell Types, Cell Neighborhoods  
① Coupling  
② Functional State  
③ Communication

**Authors**  
Christian M. Schürch, Salil S. Bhate, Graham L. Barlow, ..., Nikolay Samusik, Yury Goltsev, Garry P. Nolan

**Correspondence**  
christian.m.schuerch@gmail.com (C.M.S.), gnolan@stanford.edu (G.P.N.)

**In Brief**  
A multiplexed tissue imaging and computational analysis framework applied to colorectal cancer allows interrogation of how spatial organization of the immune tumor microenvironment is linked to clinical outcomes.

**Highlights**

- FFPE-CODEX multiplexed tissue imaging of 56 markers in 140 tissues of 35 CRC patients
- Cellular neighborhoods reveal spatial organization of the tumor microenvironment

**nature COMMUNICATIONS**

ARTICLE

<https://doi.org/10.1038/s41467-020-20803-7> OPEN

### An N-Cadherin 2 expressing epithelial cell subpopulation predicts response to surgery, chemotherapy and immunotherapy in bladder cancer

Kenneth H. Gouin<sup>1,2,3,9</sup>, Nathan Ing<sup>12,8</sup>, Jasmine T. Plummer<sup>1,2</sup>, Charles J. Rosser<sup>3,4</sup>, Bassam Ben Cheikh<sup>1,2</sup>, Catherine Oh<sup>12</sup>, Stephanie S. Chen<sup>12</sup>, Keith Syson Chan<sup>4,5</sup>, Hideki Furuya<sup>3,4</sup>, Warren G. Tourtellotte<sup>4,5,6,7</sup>, Simon R. V. Knott<sup>1,2,4,8</sup> & Dan Theodororescu<sup>3,4,5,8</sup>

**Cell Systems** Article

### Tissue schematics map the specialization of immune tissue motifs and their appropriation by tumors

Graphical abstract

**Tissue Schematics**  
Cellular Neighborhoods (CNs): regions with characteristic local processes

**Spatial Contexts**  
Fluctuations in cell states indicate local CN interactions

**Assembly Rules**  
Align tissue types, Evaluate mutations  
Good survival, Poor survival

**Authors**  
Salil S. Bhate, Graham L. Barlow, Christian M. Schürch, Garry P. Nolan

**Correspondence**  
gnolan@stanford.edu

**In brief**  
We present Tissue Schematics, a computational and conceptual framework for using high-parameter imaging data to map how tissues are formed by the assembly of cellular neighborhoods into tissue motifs. Applying our framework to human lymphoid tissues and the colorectal cancer immune-tumor microenvironment suggests mechanistic insights into the specific biology of these tissues, as well as general principles of tissue organization.

**nature COMMUNICATIONS**

ARTICLE

<https://doi.org/10.1038/s41467-020-20194-4> OPEN

### Immune cell topography predicts response to PD-1 blockade in cutaneous T cell lymphoma

Darci Phillips<sup>1,2,3,9</sup>, Magdalena Matusiak<sup>3,9</sup>, Belén Rivero Gutierrez<sup>3</sup>, Salil S. Bhate<sup>3,3,4</sup>, Graham L. Barlow<sup>1,3</sup>, Sizun Jiang<sup>1,3,5</sup>, Janos Demeter<sup>1</sup>, Kimberly S. Smythe<sup>6</sup>, Robert H. Pierce<sup>6</sup>, Steven P. Fling<sup>6</sup>, Nirasha Ramchurren<sup>6</sup>, Martin A. Cheever<sup>6</sup>, Yury Goltsev<sup>1,3</sup>, Robert B. West<sup>3</sup>, Michael S. Khodadoust<sup>7,10</sup>, Youn H. Kim<sup>2,7,10</sup>, Christian M. Schürch<sup>1,3,8,10,11</sup> & Garry P. Nolan<sup>1,3,10,12</sup>

**Blood Cancer Journal**

ARTICLE OPEN

### Lack of intrafollicular memory CD4+ T cells is predictive of early clinical failure in newly diagnosed follicular lymphoma

Patrizia Mondello<sup>1,2,13</sup>, Angelo Fania<sup>1,12</sup>, Melissa C. Larson<sup>1</sup>, Andrew L. Feldman<sup>1</sup>, Jose C. Villalobos<sup>1</sup>, Zhi-Zhang Tang<sup>1</sup>, Ilya Galkin<sup>1</sup>, Viktor Svedkolkin<sup>1</sup>, Ekaterina Postovalova<sup>1</sup>, Alexander Bagashev<sup>1</sup>, Pavel Ovcharov<sup>1</sup>, Arina Varlamova<sup>1</sup>, Sarah Huet<sup>1,9</sup>, Bruno Teisson<sup>1</sup>, Karlijn B. McGrath<sup>1</sup>, Susan Slager<sup>1</sup>, Brian K. Link<sup>1,2</sup>, Sergei Syrbu<sup>1</sup>, Anne J. Novak<sup>1</sup>, Thomas M. Habermann<sup>1</sup>, Thomas E. Witzig<sup>1</sup>, Grzegorz S. Nowakowski<sup>1</sup>, Gilles Salles<sup>1</sup>, James R. Cocharn<sup>1</sup> and Stephen M. Ansell<sup>1,14</sup>

**CSH Cold Spring Harbor Laboratory bioRxiv**

### High Resolution Single Cell Maps Reveals Distinct Cell Organization and Function Across Different Regions of the Human Intestine

John W. Hickey<sup>1\*</sup>, Winston R. Becker<sup>2\*</sup>, Stephanie A. Nevins<sup>2\*</sup>, Aaron Horning<sup>2</sup>, Almudena Espin Perez<sup>2</sup>, Roxanne Chiu<sup>2</sup>, Derek C. Chen<sup>2</sup>, Daniel Cotter<sup>2</sup>, Edward D. Esplin<sup>2</sup>, Annika K. Weimer<sup>2</sup>, Chiara Caracciolo<sup>1</sup>, Vishal Venkataramanan<sup>1</sup>, Christian M. Schürch<sup>1,7</sup>, Sarah Black<sup>1</sup>, Maria Bribi<sup>4</sup>, Kaidi Cao<sup>5</sup>, Jure Leskovic<sup>6</sup>, Zhengyan Zhang<sup>1</sup>, Shin Lin<sup>1</sup>, Teri Longacre<sup>1</sup>, Sylvia K. Plevitis<sup>1</sup>, Ying Lin<sup>5</sup>, Garry P. Nolan<sup>1@</sup>, William J. Greenleaf<sup>6@</sup>, Michael Snyder<sup>2@</sup>

**Footnotes:**  
1 Department of Pathology, Stanford School of Medicine, Stanford, CA, 94305  
2 Department of Genetics, Stanford School of Medicine, Stanford, CA, 94305  
3 Department of Biomedical Data Science, Stanford School of Medicine, Stanford, CA, 94305  
4 Department of Medicine, University of Washington, Seattle, WA, 98195  
5 Department of Surgery, Washington University, St Louis, MO 63110  
6 Department of Computer Science, Stanford University, Stanford, CA 94305, USA  
7 Department of Pathology and Neuropathology, University Hospital and Comprehensive Cancer Center Tübingen, Tübingen, Germany

\*These authors contributed equally  
@Correspondence should be addressed to: gnolan@stanford.edu; wg@stanford.edu; mpsnyder@stanford.edu

**Research Square**

Preprints are preliminary reports that have not undergone peer review. They should not be considered conclusive, used to inform clinical practice, or referenced by the media as validated information.

### A tissue atlas of ulcerative colitis to guide TNF inhibitor therapy

Aaron Mayer, Derek Holman, Utkarsh Tandon, Salil Bhate, Sunil Bodapati, Graham Barlow, Jeff Chang, Sarah Black, Erica Crenshaw, Brigid Boland, Sarah Streett, William Sandborn, John Chang, Sanjiv Gambhir, Garry Nolan, Christian Schürch, Stephan Rogalla

**Abstract**  
Ulcerative colitis is a chronic-relapsing inflammatory disease of the large intestine with a complex, multifactorial pathogenesis. TNF inhibitors are widely used to suppress immune-mediated tissue damage in ulcerative colitis patients; however, therapy failures are common. Predicting TNF inhibitor response requires an understanding of the architectural features that underlie mucosal inflammation and those responsible for resistance. Here, we used highly multiplexed immunofluorescence to uncover the spatially resolved tissue architectures underlying disease progression and treatment response in 42 tissue regions from 34 individuals. We created a tissue atlas and performed spatial analysis to identify cell-cell contacts and cellular neighborhoods. We observed that cellular functional states depend on cellular neighborhood and that a subset of inflammatory cell types and cellular neighborhoods in ulcerative colitis patients persisted even during treatment with TNF inhibitor, indicating resistant niches. A computer vision model, with no *a priori* assumptions regarding cellular architectural features, was able to predict TNF inhibitor resistance. This spatial model significantly outperformed classification models based on single-cell data alone. Our results demonstrate the value of a spatial tissue atlas as a precision medicine tool to guide treatment of patients suffering from autoimmune diseases.

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

# CODEX, Together with Sn-Sequencing, Identifies Bio-marker for Cancer

An N-Cadherin 2 expressing epithelial cell subpopulation predicts response to surgery, chemotherapy and immunotherapy in bladder cancer

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# CODEX : True Single Cell Level Protein Expression With Spatial Information


JCI insight

## Integration of spatial and single-cell transcriptomics localizes epithelial cell-immune cross-talk in kidney injury

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

Multiplexed imaging of proteins in toto with **CODEX** validates the localization of immune cell clusters inferred by Spatial Transcriptomics

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## Merger of CITE-Seq and PhenoCycler data

CITE-Seq produces multiplexed PROTEIN and mRNA expression data

CITE-Seq workflow

PhenoCycler (CODEX) results are multiplexed PROTEIN expression data

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

Resource

**Deep Profiling of Mouse Splenic Architecture with CODEX Multiplexed Imaging**

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
ScienceAdvances

## Single-cell transcriptomic analysis of mIHC images via antigen mapping

Science Advances 05 Mar 2021:  
Vol. 7, no. 10, eabc5464  
DOI: 10.1126/sciadv.abc5464

Kiya W. Govek\*, Emma C. Troisi\*, Zhen Miao, Rachael G. Aubin, Steven Woodhouse, Pablo G. Camara†

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



TABLE 1	SPECIMEN	
Tissue Type	FFPE Oropharyngeal Squamous Cell Carcinoma	FFPE Human Tonsil
Markers	103 Antibodies	103 Antibodies
Tissue imaging area	1.2 cm x 1.4 cm	1.3 cm x 1.8 cm
Fluidics time	38 cycles   44 hours	38 cycles   50 hours
Data size	24 gigabytes	46 gigabytes
Number of cells	832,823 cells	2,253,539 cells

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## RNAScope HiPlex v2

Automated HiPlex v2 Workflow on PhenoCycler-Fusion System

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