Central Research Laboratory

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Guidance for New Users FY2022 Research Equipment Section, Central Research Laboratory

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1) Overview: Usage of Research Equipment Section, CRL

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User registration process at Research Equipment Section, CRL

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- Attend this online Zoom guideline lecture
- We will provide a unique link in the Zoom chat during the lecture
- Go to the link and download "Survey form" and "Application Form"
- Fill out both forms and send them by e-mail attachment for approval
(e-mail: hqcr1@belle.shiga-med.ac.jp)

* By annual renew application, the user registration is maintained.

Details are indicated on the CRL website:

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Why is the user registration required?

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Request in the use of CRL facilities

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<Basic Rules>

- Fill out the booking sheet of each instrument that you plan to use; Write the proposed time of use, your name and belonging, and the extension number.
- Prepare any consumables, tools, chemicals, etc. for your own use.
- Strictly follow the operating rules to use instruments and analyzers.
- Fill out the user record sheet after every use (This is important for calculating the user charge).
- Take out all things brought into the working area, including waste materials (Please leave the facility in a clean condition after use).

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Cautionary notice for booking equipment

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Follow the CRL operating rules to use equipment

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Tips for users of Research Equipment Section, CRL

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- ・ Users are expected to keep rules and regulations at CRL for obtaining accurate experimental data.
- ・ CRL staffs would greatly appreciate it if the users bring a copy of an article to CRL when it is published.

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- ・ We will give a link for download site of two forms to be submitted ("survey form" and "application form for CRL User Registration") in the Zoom chat box during this course. By sending them using email, your registration will be completed.

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SUMS Industrial Safety and Health Management System

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Rules concerning Industrial Safety and Health at the CRL Facility

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SDS (Safety Data Sheet)
化学物質や化学物質が含まれる原材料などを安全に取り扱うために必要な情報を記載したもの
SDS is a form containing data regarding the properties of a particular substance. It is intended to provide workers with procedures for handling or working with that substance in a safe manner.

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How to manage and handle harmful substances

1. 有害物質の使用の中止、有害性の少ない物質への転換
Try to use less harmful substances.
2. 有害作業の改良によって有害物質発散の防止
Prevent spreading harmful substances by improving your procedure.
3. 有害物質を取り扱う設備の密閉化や自動化
Avoid leak of harmful substances from your space and introduce automated system for handling
4. 有害な行程の隔離と遠隔操作の採用
Work as far as possible from harmful substances and use remote handling equipment.
5. 局所排気装置の設置(ドラフトチャンバー)
Set up a local exhaust ventilation system (e.g. draft chamber).
6. 全体排気装置の設置
Set up a complete exhaust ventilation system.

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Example of industrial accidents in SUMS

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Introduction of Chemical Registration Information System

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Register information of all designated chemicals and high pressure gas cylinders in Chemical Registration Information System (CRIS FOREST)

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Update registered information in CRIS FOREST when designated chemicals and/or high pressure gas cylinders are used

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Introduction of Chemical Registration Information System (No.2)

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All work with formaldehyde must be performed in the fume hood of Morphological Research Joint Lab on the 4th floor. Formaldehyde Users should sign in a Formaldehyde Log after use.

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Transporting Liquid Nitrogen Container

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Emergency Procedures & Evacuation at CRL (in case of FIRE)

1. 助けを求める
Ask help without firefighting operations alone.
2. 防災監視室(内線2773)もしくはセンター職員に通報・連絡する(内線2300)
Contact immediately Security Control Center (Ext.2773) and/or the CRL office (Ext.2300).
3. 可能であれば初期消火、避難誘導、危険物の搬出に努める
Help extinguish the fire if it is at an early stage, provide evacuation instructions and bring out hazardous materials (e.g. alcohol).

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Evacuate from the fire zone if you think it dangerous.

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
How to use special laboratory designed for
pathogenic virus at CRL

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
In accordance with the SUMS Biosafety Management Rules, an application form for experiment handling pathogens, etc. should be submitted to the SUMS Biosafety Committee via the Research Promotion Division.

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How to Submit an "Application for Experiment Handling Pathogens, etc." (Form 1 or 2) designated by the SUMS Biosafety Management Rules

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Application Form 1 and 2 are available for download on the SUMS website.

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
Applications and Inquiries: The Research Promotion Division (Ext.2110)

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Application Process for Using Pathogenic Virus Lab

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
Submitting a "CRL Application Form for Handling Pathogens, etc."

Once the procedures with the Biosafety Committee are completed, a "CRL Application Form for Handling Pathogens, etc." must be submitted to CRL. The principal investigator of the experiment is required to submit a report every fiscal year.

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

Location Map of Cell Engineering Lab 4

Cell Engineering Lab 4 is on the 4th floor of CRL.

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Facility and Equipment of Pathogenic Virus Lab

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
Location Map of Pathogenic Virus Lab

Pathogenic Virus Lab is on the 4th floor of CRL.

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
4)How to use experimental animals at CRL

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Central Research Laboratory
Takuya Nakase

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Need to complete following:
1. Completed the training program in animal experiment.
2. Pass the qualification exam of animal experiment.
3. Obtain the certification on animal experiment.

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
Using animal laboratories in CRL Research Equipment Section

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
Approved animal laboratories on the 1st floor in CRL Research Equipment Section

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Law relating the recombinant DNA experiments at universities, research institutes etc.

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環境中への拡散を防止して行う使用等（大学等の実験室内）
Containment measures (Type 2 Use of LMOs)
Possible uses are those with containment measures to prevent to dispersal of LMOs in the environment (Uses in laboratories etc.)

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The government will order to suspend research.

Offender shall be punished by imprisonment or a fine, or a combination of these two.

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
An example of violation 2


This university performed virus experiments without use of a certified safety cabinet

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
Containment Measures of P1 Level Laboratory


Closed doors and windows on P1 experiment 開放厳禁
Limited access 入室制限
Washed hands
Inactivated all contaminated things

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Containment Measures of P1A(experiment with Animal) Level Laboratory

遺伝子組換え生物等が逃亡その他拡散しない構造の容器に入れること
When a LMO is taken out of the laboratory in process of an experiment, the LMO shall be put in a container of the structure that prevents it from leaking or other dispersion.

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Precautions of the recombinant DNA laboratory

共通機器を初めて使用する場合、必ずセンター職員のレクチャーを受けてください。特に遠心機は危険です
When using equipment for the first time, you must ask the CRL staffs how to operate it. Especially you must do so to use centrifuges.

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the recombinant DNA laboratory in CRL Research Equipment Section

1階 1st floor	エックス線照射室 X-ray Irradiation Lab	P1A	Yamamoto(2304) Mori(2307)
	3T MRI実験室 3-Tesla MRI lab	P1A	Nakase(2305)
3階 3rd floor	遺伝子工学実験室 Genetic Engineering lab	P2	Nakase(2305)
	小動物イメージング解析室 Small Animal Imaging Analysis Lab	P1A	Yamamoto(2304) Mori(2307) Nakase(2305)
4階 4th floor	細胞工学実験室1 Cell Engineering lab 1	P1	Mori(2307) Yamamoto(2304)
	細胞工学実験室2 Cell Engineering lab 2	P1	Mori(2307) Yamamoto(2304)
	細胞工学実験室3 Cell Engineering lab 3	P1A	Yamamoto(2304) Mori(2307)
	細胞工学実験室4 Cell Engineering lab 4	P2	Mori(2307) Yamamoto(2304)
	細胞培養室 Cell culture lab	P1	Mori(2307) Yamamoto(2304)
	病原ウイルス実験室1,2,3 Pathogenic Virus Lab 1,2,3	P2	Mori(2307) Yamamoto(2304)
レーザー顕微鏡室1 Confocal Laser Scanning Microscope Lab 1	P1A	Yamamoto(2304) Mori(2307)	

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6) Procedure for using X-ray generators

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No risk of x-ray leaking from these generators

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
Required items of education and training (4.5 hrs. ~)

- ① X-ray radiography work, ② Mechanism and safe handling of X-ray or gamma-ray irradiator, ③ Effects of ionizing radiation on the human body, ④ Laws and regulations

※Those who have already been registered as an X-ray radiography worker at the hospital may be exempted from items other than ① and ②.

学内登録申請(X線照射装置) Registration as X-ray Irradiator User

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Required items of education and training (1 hr. ~)

- ① X-ray radiography work
- ② Mechanism and safe handling of X-ray or gamma-ray irradiator

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Application for Registration as a Radiation Worker in SUMS

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Registration as an X-ray radiography worker by receiving education/training and taking a health check (blood test)

使用時はルミネスバッジを着用
Wearing a Luminess badge during use

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(2) X線照射装置
X-ray Irradiator

X線作業員
X-ray worker

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Registration as an X-ray worker by receiving education/training

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Application for Registration as a Radiation Worker in SUMS

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An application for registration of an ionizing radiation worker must be submitted to the Research Promotion Division ("Kenkyu Suisin- Ka") on the 1st floor of the D hospital bldg.

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Application for Registration as a Radiation Worker in SUMS

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Circle "2" if you wish to use an Animal X-ray Fluoroscope, and circle "3" if you wish to use an X-ray Irradiator.

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利用申込書の提出
Fill in a form for using X-ray generator

利用の承認
Acceptance for your use by CRL

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Animal X-ray Fluoroscopy room is a controlled area protected with a security system using ID card verifier.

Those wishing to start using this room need to register their ID card information prior to entry to the room for the first time.

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7) Equipment introduction #1 (Instruments on the 1st and 4th floors)

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Laboratories and instruments on the 1st floor

ヒト由来サンプルを保管する実験室、エックス線を使用した研究を行う実験室、NMR室、3T MRI実験室、工作室、ウェットラボを配置
Provided with a room where human-derived materials can be stored, laboratories for research using X-ray generators, a nuclear magnetic resonance device, a 3-Tesla MRI scanner, a mechanical workroom and Wetlab, the first floor is used for various purposes.

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Human Sample Resource Room (HSR)

- 大学内における実験・研究等に使用するヒト由来サンプルを一括して保管する施設
HSR is a facility where human-derived materials are securely and appropriately stored for research activities within the University.
- 冷凍庫の温度は温度管理システムにより監視
HSR freezers are connected to temperature alarm systems.
- 超低温槽の異常時には利用者にメール
HSR registered users can get their freezers' temperature information from any area within the University. Also, they can be contacted by these alarm systems in the event the freezer is having a temperature critical issue.

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Laboratories and instruments on the 1st floors

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NMR装置
Nuclear Magnetic
Resonance Device

工作室
Mechanical Workroom

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Laboratories and instruments on the 4th floor

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
冷却CCDカメラ付倒立型蛍光・可視顕微鏡
Inverted Fluorescence Microscope
with Cooled CCD Camera

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Microscopy Sample Preparation Equipment

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Basic rules

- 予約が必要な機器は、予め予約表に記入する
A reservation is required for the use of the equipment at CRL.
Fill out the booking sheet of each instrument that you plan to use.
- 初めて利用する時は、必ず、利用方法、操作方法、注意事項を機器担当者に問い合わせる
When using equipment at CRL for the first time, any users should ask the CRL staffs how to operate it.
- 実験室に持ち込んだものは実験終了後持ち帰る
Do not leave your items in the working area, including waste materials (Keep the facility clean after use).

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8) Equipment introduction #2 (Instruments on the 2nd and 3rd floor)

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Elix water can be used for general use for research experiments.

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Compatible to 6 to 384 well plates for measuring absorption, fluorescence, and luminescence.

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Can measure DNA, RNA, and protein from 1-2 μ L of sample in seconds.

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For protein detection on Western blot membranes.

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For quantitative measurement of gene expression.

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Gel Imaging Device FAS5 (Nippon Genetics)

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For identification and quantification of up to 100 different molecules simultaneously in a single assay.

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Focused-ultrasonicator Covaris S220 (M&S)

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For tissue homogenization, cell lysis, DNA or chromatin shearing, and other sample preparations.

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Microchip type electrophoresis device Bioanalyzer2100(Agilent)

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To analyze concentration, size, and quality of DNA, RNA, and protein.

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Instruments on the 2nd floor

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Ultracentrifuge

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High Speed Refrigerated Centrifuge

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For detection of in vivo luminescence and fluorescence in small animals.

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Instruments on the 3rd floor

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Cold chambers for biochemical experiments and storage

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
For sample storage and biochemical experiments such as chromatography at 4 °C.

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For disruption of chromatin, cells, bacteria, or tissue.

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For homogenization of tissues uniformly in a short time.

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For drying samples without denaturing under freezing temperature.

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For acquisition of anatomical, functional, physiological, and molecular level data *in vivo* in real time.

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From April 2019, a regular USB memory device is NOT allowed to transfer the data in all computers in the CRL.
We recommend to use the USB memory device with encryption function as below.

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We recommend Kingston Technology's
DataTraveler2000.

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You have 2 choices, 1) prepare DNA sequencing reaction and submit the sample for running, 2) submit a plasmid for DNA sequencing (new).

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場所：4F 405号室
Location: 4F Room405

作製 Preparation

- 光学顕微鏡用試料作製
Sample Preparations of Optical Microscope

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For more information, please check our Website, or ask us.

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Do not hesitate to ask questions when you want to use our service.

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Learn from the past accidents happened in the CRL

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The rotor came off from the spindle and rolled around in the chamber. The rotor and spindle were severely damaged.

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Two pins on the spindle were severely damaged.

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
The driving unit was damaged, while the titanium rotor installed to the spindle was intact.

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The aluminum screws of the tube caps were broken down by the centrifugal force at 50,000 rpm.
Compatible tube caps must be used for ultracentrifugation.

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Case 3: the corroded bucket of Swinging-bucket rotor was broken down during the run so that the rotor lost its balance.

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The bucket came off from the running rotor at 40,000 rpm and crashed onto the chamber wall at a velocity of over 1,600 km/h. Relative Centrifugal Field is 193,000 g.

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How to operate centrifuges safely without accidents ?

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Tips for safety centrifuge operation

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
Tips for safety operation of ultracentrifuge

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
- ▶ ローターの最高許容回転数を厳守すること。
Check the maximum allowable rotor speed of your rotor. Rotors have their operating life and, in some rotors, the maximum rotor speed should be reduced.
- ▶ 比重が1.2を超える試料を遠心するときには最高回転数を次式によって減じること。
The maximum rotor speed should be reduced using the following formula when centrifuging a solution with a density greater than 1.2.


$$\text{reduced runspeed} = \text{maximum rated speed} \times \sqrt{\frac{1.2}{\text{actual density of the tube contents}}}$$

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- ▶ 使用方法、注意事項等の説明を受けること。
User should receive the instruction of centrifuges from CRL staffs when using first time.
- ▶ 予約表への記入
Centrifuges in CRL are used on reservation basis. Please make your reservation on the sheet placed near the centrifuge.
- ▶ 利用ノートへの記入
Fill out the User-notebook after use.
- ▶ あとかたづけ
Please keep the centrifuge and the surrounding areas clean after use.

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