Central Research Laboratory

2022 Guidance for New Users Research Equipment Section Central Research Laboratory

Central Research Laboratory
Shiga University of Medical Science

2022 Guidance for New Users

Research Equipment Section, Central Research Laboratory (CRL)

1)	Overview: Usage of Research Equipment Section, CRL	1
2)	Industrial safety and health	4
3)	How to use the pathogenic virus laboratory	7
4)	How to use experimental animals at CRL	10
5)	How to use the recombinant DNA laboratory	12
6)	Procedures for using X-ray generators	14
7)	Equipment introduction #1 (Instruments on the 1st and 4th floors)	17
8)	Equipment introduction #2 (Instruments on the 2nd and 3rd floors)	20
9)	Special service rendered by CRL	24
10)	Safety operation for centrifuges	26

Appendix: Equipment List / CRL Floor Map / Staff List

***** How to Register *****

After taking the guidance, you can download the forms to be submitted ("survey form" and "application form for CRL User Registration"). By sending them using email, your registration will be completed.

Once the registration has been processed...

Registered users can continue their registration by simply going through the renewal procedure (submitting the renewal application form) from the following year.

Central Research Laboratory, SUMS
Tel: 077-548-2300
E-mail: hqcrl@belle.shiga-med.ac.jp
URL: <http://wwwcrl.shiga-med.ac.jp/index.html>

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Guidance for New Users FY2022 Research Equipment Section, Central Research Laboratory

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1) Overview: Usage of Research Equipment Section, CRL

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User registration process at Research Equipment Section, CRL

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1. Attend this online Zoom guideline lecture
2. We will provide a unique link in the Zoom chat during the lecture
3. Go to the link and download "Survey form" and "Application Form"
4. Fill out both forms and send them by e-mail attachment for approval
(e-mail: hqcr1@belle.shiga-med.ac.jp)

* By annual renew application, the user registration is maintained.

Details are indicated on the CRL website:

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Why is the user registration required?

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Request in the use of CRL facilities

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<Basic Rules>

1. Fill out the booking sheet of each instrument that you plan to use; Write the proposed time of use, your name and belonging, and the extension number.
2. Prepare any consumables, tools, chemicals, etc. for your own use.
3. Strictly follow the operating rules to use instruments and analyzers.
4. Fill out the user record sheet after every use (This is important for calculating the user charge).
5. Take out all things brought into the working area, including waste materials (Please leave the facility in a clean condition after use).

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Cautionary notice for booking equipment

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1. Write the extension number on the booking sheet.
Users should communicate with other users and the CRL staffs.
2. Do not make booking for inappropriately long period of time.
The equipment booking may be canceled if the equipment is not used after 30 minutes of proposed starting time, especially for equipment that is commonly used by many users.
3. Follow the booking rules for each equipment.
4. Follow the CRL rules and strictly observe the start/end time for use.

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Follow the CRL operating rules to use equipment

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- ・When using equipment at CRL for the first time, **any users MUST ask the CRL staffs** how to operate it.
- ・Online operation manuals are available on the CRL website (http://www.crl.shiga-med.ac.jp/home/kiki_bumon/g_book/contents.html). If necessary, print them out for your reference to use the equipment.

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Tips for users of Research Equipment Section, CRL

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-
- ・ Users are expected to keep rules and regulations at CRL for obtaining accurate experimental data.
 - ・ CRL staffs would greatly appreciate it if the users bring a copy of an article to CRL when it is published.

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-
- ・ We will give a link for download site of two forms to be submitted ("survey form" and "application form for CRL User Registration") in the Zoom chat box during this course. By sending them using email, your registration will be completed.

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SDS (Safety Data Sheet)
化学物質や化学物質が含まれる原材料などを安全に取り扱うために必要な情報を記載したもの
SDS is a form containing data regarding the properties of a particular substance. It is intended to provide workers with procedures for handling or working with that substance in a safe manner.

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Introduction of Chemical Registration Information System

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Register information of all designated chemicals and high pressure gas cylinders in Chemical Registration Information System (CRIS FOREST)

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Update registered information in CRIS FOREST when designated chemicals and/or high pressure gas cylinders are used

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Introduction of Chemical Registration Information System (No.2)

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All work with formaldehyde must be performed in the fume hood of Morphological Research Joint Lab on the 4th floor. Formaldehyde Users should sign in a Formaldehyde Log after use.

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Transporting Liquid Nitrogen Container

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Emergency Procedures & Evacuation at CRL (in case of FIRE)

1. 助けを求める
Ask help without firefighting operations alone.
2. 防災監視室(内線2773)もしくはセンター職員に通報・連絡する(内線2300)
Contact immediately Security Control Center (Ext.2773) and/or the CRL office (Ext.2300).
3. 可能であれば初期消火、避難誘導、危険物の搬出に努める
Help extinguish the fire if it is at an early stage, provide evacuation instructions and bring out hazardous materials (e.g. alcohol).

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
Evacuate from the fire zone if you think it dangerous.

Emergency Procedures & Evacuation at CRL (in case of ACCIDENT)

1. **まず自分の安全を確保し二次災害を防ぐ**
Make sure you are safe and then prevent the secondary accident.
2. **怪我人がいる場合は適切な処置を行う**
Help injured persons if it is possible.
3. **センター職員に連絡する(内線2300)**
Notify the condition to the CRL staffs (Ext.2300) immediately.
4. **後日、事故報告書を支援センターに提出**
Submit an accident report to CRL.

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
How to use special laboratory designed for pathogenic virus at CRL

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
In accordance with the SUMS Biosafety Management Rules, an application form for experiment handling pathogens, etc. should be submitted to the SUMS Biosafety Committee via the Research Promotion Division.

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How to Submit an "Application for Experiment Handling Pathogens, etc."
(Form 1 or 2) designated by the SUMS Biosafety Management Rules

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Application Form 1 and 2 are available for download on the SUMS website.

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
Applications and Inquiries: The Research Promotion Division (Ext.2110)

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Application Process for Using Pathogenic Virus Lab

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
Submitting a "CRL Application Form for Handling Pathogens, etc."

Once the procedures with the Biosafety Committee have been completed, a "CRL Application Form for Handling Pathogens, etc." must be submitted to CRL. The principal investigator of the experiment is required to submit a report every fiscal year.

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

Location Map of Cell Engineering Lab 4

Cell Engineering Lab 4 is on the 4th floor of CRL.

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Facility and Equipment of Pathogenic Virus Lab

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
Location Map of Pathogenic Virus Lab

Pathogenic Virus Lab is on the 4th floor of CRL.

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Layout of Pathogenic Virus Lab

Adenovirus
Adenovirus
Influenza virus
Influenza virus

Biosafety cabinet and CO2 incubator for Adenovirus and Influenza virus

Retrovirus
Retrovirus
Lentivirus
Lentivirus

Biosafety cabinet and CO2 incubator for Retrovirus and Lentivirus

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
4)How to use experimental animals at CRL

Central Research Laboratory
Takuya Nakase

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- Need to complete following:
1. Completed the training program in animal experiment.
 2. Pass the qualification exam of animal experiment.
 3. Obtain the certification on animal experiment.

修正
Revise

承認
Approve

Obtain approval for registration to
CRL Research Equipment Section

Obtain approval for registration to
CRL Research Equipment Section

Using animal laboratories in CRL Research Equipment Section

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
Approved animal laboratories on the 1st floor in CRL Research Equipment Section

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Law relating the recombinant DNA experiments at universities, research institutes etc.

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環境中への拡散を防止して行う使用等（大学等の実験室内）
Containment measures (Type 2 Use of LMOs)
Possible uses are those with containment measures to prevent to dispersal of LMOs in the environment (Uses in laboratories etc).

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The government will order to suspend research.

Offender shall be punished by imprisonment or a fine, or a combination of these two.

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
An example of violation 2


This university performed virus experiments without use of a certified safety cabinet

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
P1レベル拡散防止措置 (抜粋)
Containment Measures of P1 Level Laboratory


Closed doors and windows on P1 experiment

Limited access

Washed hands

Inactivated all contaminated things

P1レベル拡散防止措置 (抜粋)
Containment Measures of P1A(experiment with Animal) Level Laboratory
遺伝子組換え生物等が逃亡その他拡散しない構造の容器に入れること
When a LMO is taken out of the laboratory in process of an experiment, the LMO shall be put in a container of the structure that prevents it from leaking or other dispersion.
遺伝子組換え生物等の種類ごとに識別することができる措置を講ずること
Take a measure that is capable of identifying modified animals.

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Precautions of the recombinant DNA laboratory

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6) Procedure for using X-ray generators

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
No risk of x-ray leaking from these generators

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

利用登録の手順 Registration procedure

Registration of Animal X-ray Fluoroscope User

Taking education and training → Registration as an Animal X-ray Fluoroscope User (X-ray radiography worker) → Blood test

※Those who have already been registered as an X-ray radiography worker at the hospital may be exempted from items other than ① and ②.

Registration as X-ray Irradiator User

Taking education and training → Registration as an X-ray Irradiator User (X-ray worker)

- ① X-ray radiography work
- ② Mechanism and safe handling of X-ray or gamma-ray irradiator

Application for Registration as a Radiation Worker in SUMS

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使用前に教育訓練
Registration as an X-ray
worker by receiving
education/training

Application for Registration as a Radiation Worker in SUMS

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Application for Registration as a Radiation Worker in SUMS

[illegible]

An application for registration of an ionizing radiation worker must be submitted to the Research Promotion Division ("Kenkyu Suisin- Ka") on the 1st floor of the D hospital bldg.

Application for Registration as a Radiation Worker in SUMS

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Circle "2" if you wish to use an Animal X-ray Fluoroscope, and circle "3" if you wish to use an X-ray Irradiator.

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利用申込書の提出
Fill in a form for using X-ray generator

利用の承認
Acceptance for your use by CRL

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Animal X-ray Fluoroscopy room is a controlled area protected with a security system using ID card verifier.

Those wishing to start using this room need to register their ID card information prior to entry to the room for the first time.

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Please contact the CRL staff for more information.

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7) Equipment introduction #1 (Instruments on the 1st and 4th floors)

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Laboratories and instruments on the 1st floor

ヒト由来サンプルを保管する実験室、エックス線を使用した研究を行う実験室、NMR室、3T MRI実験室、工作室、ウェットラボを配置
Provided with a room where human-derived materials can be stored, laboratories for research using X-ray generators, a nuclear magnetic resonance device, a 3-Tesla MRI scanner, a mechanical workroom and Wetlab, the first floor is used for various purposes.

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Human Sample Resource Room (HSR)

- 大学内における実験・研究等に使用するヒト由来サンプルを一括して保管する施設
HSR is a facility where human-derived materials are securely and appropriately stored for research activities within the University.
- 冷凍庫の温度は温度管理システムにより監視
HSR freezers are connected to temperature alarm systems.
- 超低温槽の異常時には利用者にメール
HSR registered users can get their freezers' temperature information from any area within the University. Also, they can be contacted by these alarm systems in the event the freezer is having a temperature critical issue.

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
Laboratories and instruments on the 1st floors

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NMR装置
Nuclear Magnetic
Resonance Device

工作室
Mechanical Workroom

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Laboratories and instruments on the 4th floor

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


冷却CCDカメラ付倒立型蛍光・可視顕微鏡
Inverted Fluorescence Microscope
with Cooled CCD Camera

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Microscopy Sample Preparation Equipment

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Basic rules

- 予約が必要な機器は、予め予約表に記入する
A reservation is required for the use of the equipment at CRL.
Fill out the booking sheet of each instrument that you plan to use.
- 初めて利用する時は、必ず、利用方法、操作方法、注意事項を機器担当者に問い合わせる
When using equipment at CRL for the first time, any users should ask the CRL staffs how to operate it.
- 実験室に持ち込んだものは実験終了後持ち帰る
Do not leave your items in the working area, including waste materials (Keep the facility clean after use).

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8) Equipment introduction #2 (Instruments on the 2nd and 3rd floor)

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Compatible to 6 to 384 well plates for measuring absorption, fluorescence, and luminescence.

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Can measure DNA, RNA, and protein from 1-2 μ L of sample in seconds.

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

For protein detection on Western blot membranes.

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For quantitative measurement of gene expression.

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Gel Imaging Device FAS5 (Nippon Genetics)

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
For identification and quantification of up to 100 different molecules simultaneously in a single assay.

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Focused-ultrasonicator Covaris S220 (M&S)

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For tissue homogenization, cell lysis, DNA or chromatin shearing, and other sample preparations.

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Microchip type electrophoresis device Bioanalyzer2100(Agilent)

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To analyze concentration, size, and quality of DNA, RNA, and protein.

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Instruments on the 2nd floor

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Ultracentrifuge

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High Speed Refrigerated Centrifuge

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For detection of in vivo luminescence and fluorescence in small animals.

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Instruments on the 3rd floor

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Cold chambers for biochemical experiments and storage

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
For sample storage and biochemical experiments such as chromatography at 4 °C.

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For disruption of chromatin, cells, bacteria, or tissue.

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For homogenization of tissues uniformly in a short time.

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For drying samples without denaturing under freezing temperature.

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High-resolution Ultrasonic Echo Imaging System for Small Animals

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For acquisition of anatomical, functional, physiological, and molecular level data *in vivo* in real time.

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From April 2019, a regular USB memory device is NOT allowed to transfer the data in all computers in the CRL.

We recommend to use the USB memory device with encryption function as below.

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We recommend Kingston Technology's
DataTraveler2000.

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You have 2 choices, 1) prepare DNA sequencing reaction and submit the sample for running, 2) submit a plasmid for DNA sequencing (new).

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
場所：4F 405号室
Location: 4F Room405

作製 Preparation


- 光学顕微鏡用試料作製
Sample Preparations of Optical Microscope

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For more information, please check our Website, or ask us.

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Do not hesitate to ask questions when you want to use our service.

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Learn from the past accidents happened in the CRL

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
The rotor came off from the spindle and rolled around in the chamber. The rotor and spindle were severely damaged.

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
Two pins on the spindle were severely damaged.

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
The driving unit was damaged, while the titanium rotor installed to the spindle was intact.

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The aluminum screws of the tube caps were broken down by the centrifugal force at 50,000 rpm.
Compatible tube caps must be used for ultracentrifugation.

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
Case 3: the corroded bucket of Swinging-bucket rotor was broken down during the run so that the rotor lost its balance.

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The bucket came off from the running rotor at 40,000 rpm and crashed onto the chamber wall at a velocity of over 1,600 km/h.
Relative Centrifugal Field is 193,000 g.

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How to operate centrifuges safely
without accidents ?

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Tips for safety centrifuge operation

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

Tips for safety operation of ultracentrifuge

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

- ▶ 使用方法、注意事項等の説明を受けること。
User should receive the instruction of centrifuges from CRL staffs when using first time.
- ▶ 予約表への記入
Centrifuges in CRL are used on reservation basis. Please make your reservation on the sheet placed near the centrifuge.
- ▶ 利用ノートへの記入
Fill out the User-notebook after use.
- ▶ あとかたづけ
Please keep the centrifuge and the surrounding areas clean after use.

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Equipment List (Research Equipment Section, Central Research Laboratory)

BIOCHEMISTRY	EQUIPMENT NAME	FEE (JPY)	LOCATION
Microplate Reader	TECAN/Infinite M200・M PLEX	100/use	2F Biochemical Analysis Lab 1
Micro Mass Spectrometer	Thermo Fisher/NanoDrop One, NanoDrop ND-1000	—	2F Biochemical Analysis Lab 1
Recording Spectrophotometer	SHIMADZU/UV-3100PC	—	2F Protein Engineering Lab
Fourier Transform Infrared Spectrometer	SHIMADZU/IRPrestige-21	—	2F Protein Engineering Lab
Fluorescence Spectrometer	HITACHI/F-4500	—	2F Protein Engineering Lab
Circular Dichroism spectrometer	JASCO/J-600	—	2F Biochemical Analysis Lab 1
Gel Imaging System	NIPPON Genetics/FAS-IV・FAS5	—	2F Molecular Biology Analysis Lab
Nuclear Magnetic Resonance Spectroscopy	JEOL/JNM-ECZ400S	900/h (Day Max 3,000)	1F NMR Lab
High Performance Liquid Chromatograph	SHIMADZU/HPLC system	200/use	2F Biochemical Analysis Lab 1
HPLC AKTA System	GE	200/use	2F Protein Engineering Lab
Preparative Ultracentrifuge	Beckman/Optima L-90K	100/10 min (5,000/More than 500min)	2F Centrifuge Room
Tabletop Ultracentrifuge	Beckman/Optima TLX	100/10 min (5,000/More than 500min)	2F Centrifuge Room
High Speed Refrigerated Centrifuge	Beckman/HP-25・HP-26	200/use	2F Centrifuge Room
Hybrid High Speed Refrigerated Centrifuge	KUBOTA/6200	—	2F Centrifuge Room
Centrifugal Concentrator/Speedvac Concentrator	Thermo Fisher/SPD1010	—	2F Centrifuge Room
High Speed Refrigerated Centrifuge (P2 Lab)	Beckman/Avanti J-E	200/use	2F Gene Engineering Lab P2
High Speed Refrigerated Microcentrifuge	TOMY/MX-305	—	2F Gene Engineering Lab P2
	HITACHI/CF16RN	—	3F Biological Sample Analysis Lab
Plate Centrifuge	KUBOTA/PlateSpinII	—	2F Molecular Biology Analysis Lab
Hybrid Centrifuge	Beckman/GS-15R	—	2F Gene Engineering Lab P2
Gene Transfer System	Bio-Rad/Gene Pulser II	—	2F Gene Engineering Lab P2
Deep Freezer	NIHON FREEZER/CLN-52U	6,000/y	2F Gene Engineering Lab P2
Shaking Incubator	Sanki Seiki/SCS-R・SCS-12R, TAITEC/BR-30L, BR-30LF, BR-	—	2F Gene Engineering Lab P2
Sample Heating and Stirring Equipment	GL Sciences/Reacti-Therm 1030-44010	—	3F Analysis Room
Vacuum Freeze Drying Equipment	LABCONCO/LL-1	200/d	3F Cell homogenization Room
Ultrasonic Homogenizer	BRANSON/Sonifier Model 450-Advanced	—	3F Cell homogenization Room
Cap Horn Homogenizer	MISONIX/ASTRASON MODEL XL2020	—	3F Cell homogenization Room
Polytron Homogenizer	KINEMATICA/MODEL PT1200E	—	3F Cell homogenization Room
Chromatochamber	NIHON FREEZER/MC-30EF3	500/rack/m	2F Frozen Storage Room
Gas Chromatograph (EI)	SHIMADZU/GCMS-QP5050A	1,000/use	2F MS Spectrometry Lab
Inductively Coupled Plasma Mass Spectrometer	Agilent/7700e ICP-MS	Basic charge 800 +1,400/h	2F MS Spectrometry Lab

BIOCHEMISTRY	EQUIPMENT NAME	FEE (JPY)	LOCATION
X-ray Diffractometer	Philips/MPD-1880	1,000/use	1F X-ray Irradiation Lab
DNA Sequencer	ABI/PRISM 3130XL GeneticAnalyzer	200/sample	3F Gene Analysis Lab
Acoustic Solubilizer	Covaris/S220	100/use	2F Molecular Biology Analysis Lab
Real Time Quantitative PCR (LightCycler)	Roche/LightCycler480 System II	100/use	2F Molecular Biology Analysis Lab
Real Time Quantitative PCR (Plate)		7,600/set	2F Molecular Biology Analysis Lab
Thermal Cycler	ABI/GeneAmp PCR System 2720, Roche/LightCycler	—	2F Molecular Biology Analysis Lab
	Thermo Fisher/MiniAmp Plus	—	3F Molecular Biology Analysis Lab
	Applied Biosystems/GeneAmp PCR System 9600	—	3F Analysis Room
Lumino Image Analyzer	FUJIFILM/LAS-4000	100/use	2F Molecular Biology Analysis Lab
Chemiluminescence Imaging System	Vilber-Lourmat/FUSION	100/use	2F Molecular Biology Analysis Lab
Electrophoresis Apparatus	Pharmacia/IPG Phor	500/use	2F MS Spectrometry Lab
Bioanalyzer	Agilent Technologies/Agilent 2100 Bioanalyzer Limited	100/use	2F Molecular Biology Analysis Lab
Suspension Array System	Bio-Rad/Bio-Plex200	3,000/use	2F Biochemical Analysis Lab 1
MOE (Molecular Operating Environment)	Chemical Computing Group/MOE	—	3F Analysis Room

CELL CULTIVATION	EQUIPMENT NAME	FEE (JPY)	LOCATION
Electroporator	LONZA/4D-Nucleofector, amaxa/Nucleofector	—	4F Cell Culture Lab
CO2 Incubator	YAMATO/IP400, SANYO/MCD-175, ESPEC/BNA-111, HIRASAWA/CPD- 2701	45,000/y + CO2 cost	4F Cell Culture Lab
CO2		6,000/cylinder	4F Cell Culture Lab
Clean Bench (Laminar Flow Cabinet)	ESCO/Airstream PRO ESC-BCB-4A7, Showa Kagaku/S-1801WBV, etc.	—	4F Cell Culture Lab
Flexercell	FLEXERCELL/FLEX-I	2,000/use	4F Cell Culture Lab
Microinjection System	Eppendorf/InjectMan NI2, OLYMPUS/IX-73	100/use	4F Cell Engineering Lab 3
Cell Analyzer	BD/FACSCalibur・FACSCanto II	1,500/use	4F Cell Engineering Lab 2
	BD/LSRFortessa	1,500/use	4F Cell Engineering Lab 4
Cell Sorter	BD/FACSAria	5,000/use (Analysis)1,500/use	4F Cell Engineering Lab 1
	BD/FACSAria Fusion	5,000/use (Analysis)1,500/use	4F Cell Engineering Lab 4
Cell Dynamic Analysis System	GE/EZ-TAXIScan	unfixed	4F Cell Culture Lab
Extracellular Flux Analyzer	Seahorse/XFe24	100/use	4F Cell Culture Lab
Magnetic Cell Separator	Miltenyi/autoMACS Pro Separator	unfixed	4F Cell Culture Lab
Deep Freezer (Frozen Strage Room) Container Basket	REVCO/ULT-1490-3J-D30, etc.	6,000/y	3F Frozen Strage Room
Deep Freezer (Frozen Strage Room) Drawer Case (Large)	REVCO/ULT-1490-3J-D30, etc.	4,000/y	3F Frozen Strage Room
Large Liquid Nitrogen Sample Storage Container	CHART/MVE Heco1536P- 190AF・GB	550/m	4F Workroom
Biosafety Cabinet	ESCO/ClassII Type A2	—	4F Pathogenic Virus Lab

CELL CULTIVATION	EQUIPMENT NAME	FEE (JPY)	LOCATION
Biosafety Cabinet	Thermo/1387	—	4F Pathogenic Virus Lab
Preparative Ultracentrifuge	Beckman/Optima L-90K	100/10 min (5,000/More than 500min)	4F Pathogenic Virus Lab
High Speed Refrigerated Centrifuge	Beckman/Avanti J-E	200/use	4F Pathogenic Virus Lab
Tabletop Centrifuge	KUBOTA/KS-5000	—	4F Pathogenic Virus Lab
Tabletop Refrigerated Centrifuge	Beckman/Allegrax-30R	—	4F Cell Culture Lab
Ethylene Oxide Gas Sterilizer	ELK CORPORATION/ SA-360ECO	3,000/use	4F Frozen Strage Room
X-ray Irradiator	HITACHI/MBR-1520R	1,000/use	1F X-ray Irradiation Lab
	AcroBio/RS320	1,000/use	Research Center for Animal Life Science


MORPHOLOGY	EQUIPMENT NAME	FEE (JPY)	LOCATION
Transmission Electron Microscope	HITACHI/H-7500	500/h + film cost (Market Price)	4F Electron Microscope Lab
Scanning Electron Microscope	JEOL/JSM-7505FA	1,000/30min	4F Electron Microscope Lab
Ultramicrotome	Reichert/ULTRACUTE	500/use	4F Microscopy Sample Preparation Room 1
Cooled CCD Camera	NIPPON ROPER/CoolSNAP-HQ	300/use	4F Electron Microscope Lab
Color CCD Camera	NIPPON ROPER/Retiga2000R	100/use	4F Cell Culture Room 4F Image Processing Room
High-resolution Fluorescence Microscope	GE/DeltaVision Elite	200/use	4F Microscope Lab 1
All-in-One Fluorescence Microscope	KEYENCE/BIOREVO BZ-9000	—	4F Cell Culture Lab
Inverted Fluorescence Microscope	OLYMPUS/IX83	100/use	4F Microscopy Sample Preparation Room 2
Confocal Laser Scanning Microscope	ZEISS/LSM-510 META	1,000/use	4F Microscopy Sample Preparation Room 2
	NIKON/C1si	1,000/use	4F Confocal Laser Scanning Microscope Lab 2
	OLYMPUS/FV1000-D	1,000/use	4F Microscope Lab 2
	LEICA/TCS SP8 X	1,000/use	4F Confocal Laser Scanning Microscope Lab 3
Live Cell Imaging Microscope	Andor/Dragonfly200	300/h (Day Max 5,000)	4F Confocal Laser Scanning Microscope Lab 1
Microtome	LEICA/SM2010R	100/use	4F Microscopy Sample Preparation Room 2
Tissue Processor	SAKURA/Tissue-Tek VIP 6-J0	1,000/use	4F Microscopy Sample Preparation Room 2
Tissue Embedding Console	SAKURA/Tissue-Tek 6	100/piece	4F Microscopy Sample Preparation Room 2
Cryostat	LEICA/CM3050 S	500/use	4F Microscopy Sample Preparation Room 2
Sample Preparation	Paraffin embedding	500/use	4F Microscopy Sample Preparation Room 1
	Sectioning	300/slide	4F Microscopy Sample Preparation Room 1
	H&E Staining	300/slide	4F Microscopy Sample Preparation Room 1
3D/4D Image Analysis and Visualization Software	BITPLANE/Imaris	100/use	4F Image Processing Room

PHYSIOLOGY / ANIMAL	EQUIPMENT NAME	FEE (JPY)	LOCATION
Microdialysis	EICOM/HTEC-500	1,000/use	2F Biochemical Analysis Lab 1
Animal Blood Gas Analyzer	ABAXIS/VetScan VS2	—	3F Biological Sample Analysis Lab
Automatic Blood Cell Counter	NIHON KOHDEN/ Celltac α MEK-6550	200/sample	3F Biological Sample Analysis Lab
Oxidation Redox Analyzer	WISMERLL/REDOXLIBRA	100/use	3F Biological Sample Analysis Lab
Laboratory Animal X-ray Fluoroscope	TOSHIBA Medical/Plessart50	1,000/use	1F Animal X-ray Fluoroscopy Lab
Langendorff System	ADInstruments	1,000/use	4F Cell Engineering Lab 3
High-resolution Ultrasonic Echo Imaging System for Small Animals	VisualSonics/Vevo2100	1,500/use	3F Small Animal Imaging Analysis Lab
In vivo Imaging System for Small Animals	Vilber-Lourmat/NEWTON 7.0	200/h	3F Small Animal Imaging Analysis Lab
3-Tesla MRI Scanner	SIEMENS Healthineers/ MAGNETOM Verio dot 3T	400,000/y 3,000/d, 1,000/h	1F 3T MRI Lab
4.7-Tesla MRI Scanner	BRUKER/BioSpec47/40USR	400,000/y 3,000/d, 1,000/h	Research Center for Animal Life Science
7-Tesla MRI Scanner	Varian/Agilent Technologies/ UNITY INOVA 7T/400/AS	240,000/y 3,000/d, 1,000/h	MR Research Facility


ADDITIONAL SERVICE	EQUIPMENT NAME	FEE (JPY)	LOCATION
Elix Water	MILLIPORE/Milli-Q Integral 5	10/L	2F Biochemical Analysis Lab 4F Workroom
Milli-Q Water	MILLIPORE/Milli-Q Integral 5	80/L	2F Biochemical Analysis Lab 4F Workroom
Ice Machine	Scotsman/F0522B Hoshizaki/FM-340AK-SA	—	2F Biochemical Analysis Lab 4F Workroom
Liquid Nitrogen	100 liter, 50 liter tank	400/L (Market Price)	4F Workroom
Oxygen Gas	purity 99.5%, 4.0Kg/cm ²	3,630/cylinder (Market Price)	1F Bombe Room 1F Wet Lab
Compressed Air		5,500/cylinder (Market Price)	1F Wet Lab
Human Sample Storing	REVCO, NIHON FREEZER	2,000/freezer/m	1F Human Sample Resource Room
Space Charge		2,000/3.3m ² /m	
Equipment Rental		200/w	

WET LAB	EQUIPMENT NAME	FEE (JPY)	LOCATION
Wet Lab Space Rental		2,000/3.3m ² /m	1F Wet Lab

1F·2F of Central Research Laboratory

3F·4F of Central Research Laboratory


Application Procedures for New User Registration

Those who wish to use the facilities and equipment of Central Research Laboratory (CRL) must be registered users. To be the user, they are required to take the Guidance for New Users managed by CRL. Depending on the nature of the experiments they plan to conduct, other application procedures may be required in addition to the user registration.

They must complete the necessary procedures according to the section they wish to use and the type of experiments they are planning to conduct.

Once you have completed the registration ···

Each piece of equipment at the CRL has its own rules on how to reserve and use it. If you use equipment for the first time, please be sure to ask the staff member in charge of the equipment how to use it before using it, rather than having your lab members teach you how to use it.

Items that Require Application to the University

- Recombinant DNA Experiments: Application to the SUMS Genetic Recombination Experimentation Safety Committee (Research Promotion Division)
- Experiments Handling pathogenic microorganisms: Application to the SUMS Biosafety Committee (Research Promotion Division)
- Experimental Animals: Application to the SUMS Animal Experimentation Committee (Research Center for Animal Life Science)
- Experiments with Animal X-ray Fluoroscope: Registration as an X-ray photographing worker (Research Promotion Division)
- Experiments with X-ray Irradiator: Registration as an X-ray worker (Research Promotion Division)
- Human Subjects Research: Application to the SUMS Ethics Committee

Please follow the instructions below when completing “Application Form for CRL User Registration”

The approval number and other necessary information must be entered on the form, and the documents such as an experiment plan must be attached, if necessary, to conduct the following experiments at CRL:

- (1) **Recombinant DNA experiments**: Enter the approval number. Persons using the Gene Engineering Lab must register their ID card information to the entry/exit managing system prior to using the Lab.
- (2) **Experiments handling pathogenic microorganisms**: Enter the approval number and submit the prescribed documents. Persons using the Pathogenic Virus Lab and/or the Cell Engineering Lab #4 must register their ID card information to the entry/exit managing system prior to using the Lab(s).
- (3) **Animal experiments**: Enter the approval number. Persons using the **Wet Lab** must submit the copy of the experiment plan submitted to the Research Center for Animal Life Science.
- (4) **Experiments using X-ray generators** (Animal X-ray Fluoroscope, X-ray Irradiator, X-ray Diffractometer): Take the X-ray training course.
- (5) **Human Sample Resource Room**: Submit the prescribed form. Persons using the HSR must register their ID card information to the entry/exit managing system prior to using the room.

【Caution】 Matters to be entered, documents to be submitted

Items	Approval number	Registration of ID card	Document to submit	Documents to submit
Recombinant DNA experiments	●	● (P2)		
Handling pathogenic microorganisms	●	●	●	①CRL Application Form for Handling Pathogens, ②Application for Handling Pathogens (Form 1 or 2)
Experimental Animals	●			
X-ray generators		●		
Wet Lab	●		●	Experiment plan for use of the Research Center for Animal Life Science
Human Sample Resource Room		●		HSR Access Application

Contact Number & E-mail Address

Research Equipment Section of Central Research Laboratory

○ Director of CRL

Job title	Name	Ext.	E-mail
Professor	OGITA, Hisakazu	2161	hogita@belle.shiga-med.ac.jp

○ Faculty of CRL

Job title	Name	Ext.	E-mail
Associate Professor	ASAHINA, Kinji	2301	asahina@belle.shiga-med.ac.jp

○ Technical Support Staff

Job title	Name	Ext.	E-mail
Technical Staff	URUSHIYAMA, Noboru	2302	urushi@belle.shiga-med.ac.jp
Technical Staff	YAMAMOTO, Takefumi	2304	tyama@belle.shiga-med.ac.jp
Technical Staff	MORI, Yasuhiro	2307	moriyasu@belle.shiga-med.ac.jp
Technical Staff	KOYAMA, Yukiko	2303	ykoyama@belle.shiga-med.ac.jp
Technical Staff	NAKASE, Takuya	2305	tnakase@belle.shiga-med.ac.jp

○ Administrative Staff

Job title	Name	Ext.	E-mail
Administrative Staff	YAMAKAWA, Nobuko	2300	nyama@belle.shiga-med.ac.jp

○ Central Research Laboratory

Location	Ext.	E-mail
2F Central Research Laboratory office	2300	hqcr1@belle.shiga-med.ac.jp

More information on CRL is available at the CRL web site (<http://wwwcrl.shiga-med.ac.jp/index.html>) where you can get a list of installed equipment, online edition of the operation manuals, etc. and download various request/application forms.

EXAMPLE

Take extra care not to make mistake or omit anything

Central Research Laboratory (CRL), Research Equipment Section

2022 New User Registration Form

Date:

I will abide by rules and regulations established by CRL.

Affiliation	
Position	
Name	
Staff (Student) ID number	
Phone number	
e-mail address	@belle.shiga-med.ac.jp
Supervisor's name	

Brief description of your research projects :

Please check boxes you plan to use

<input type="checkbox"/> Spectrophotometer	<input type="checkbox"/> Bio-Plex suspension array	<input type="checkbox"/> Live imaging microscope	<input type="checkbox"/> Heart perfusion system
<input type="checkbox"/> Mass spectrometer	<input type="checkbox"/> MOE program	<input type="checkbox"/> Fluorescence microscope	<input type="checkbox"/> Small animal imaging
<input type="checkbox"/> Nuclear magnetic resonance	<input type="checkbox"/> Chemilumi imaging	<input type="checkbox"/> CO ₂ incubator	<input type="checkbox"/> In vivo imaging
<input type="checkbox"/> Temperature control centrifuge	<input type="checkbox"/> PCR/qPCR	<input type="checkbox"/> MACS	<input type="checkbox"/> VetScan
<input type="checkbox"/> Ultracentrifuge	<input type="checkbox"/> Gel imaging	<input type="checkbox"/> Flow cytometer	<input type="checkbox"/> X ray imaging
<input type="checkbox"/> Lyophilizer	<input type="checkbox"/> Chromatography chamber	<input type="checkbox"/> FACS	<input type="checkbox"/> X ray exposure
<input type="checkbox"/> Vaccume concentrator	<input type="checkbox"/> Electronmicroscope	<input type="checkbox"/> EZ-TAXIScan	<input type="checkbox"/> MRI
<input type="checkbox"/> Micro-dialysis analyzer	<input type="checkbox"/> EM sample preparation	<input type="checkbox"/> Nucleofector	<input type="checkbox"/> Millipore water
<input type="checkbox"/> Homogenizer	<input type="checkbox"/> Light microscope	<input type="checkbox"/> Flux analyzer	<input type="checkbox"/> Liquid nitrogen
<input type="checkbox"/> DNA sequencer	<input type="checkbox"/> Microscope sample prep	<input type="checkbox"/> Deep freezer	<input type="checkbox"/> Oxygen gas
<input type="checkbox"/> Bio analyzer	<input type="checkbox"/> Image analyzer	<input type="checkbox"/> Liquid nitrogen dewar	<input type="checkbox"/> Handicraft tools
<input type="checkbox"/> Covaris solubilizer	<input type="checkbox"/> Laser microscope	<input type="checkbox"/> Autoclave	<input type="checkbox"/> Human sample storage

Please check boxes if your experiment at CRL is involved in recombinant DNA, pathogenic viruses and microorganisms, or animal use

※ New users will need to submit an additional application for entering these room by a card key

<input type="checkbox"/> Recombinant DNA (<input type="checkbox"/> New · <input type="checkbox"/> Renewal)	Approved number :		Supervisor's name : (Affiliation :)	
	Approved room	<input type="checkbox"/> Genetic Engineering Lab (※) <input type="checkbox"/> X-ray Irradiation Lab <input type="checkbox"/> 3T MRI Lab <input type="checkbox"/> Cell Culture Lab <input type="checkbox"/> Animal Imaging Lab <input type="checkbox"/> Laser Microscope Lab 1 <input type="checkbox"/> Cell Engineering Lab 1 <input type="checkbox"/> Cell Engineering Lab 2 <input type="checkbox"/> Cell Engineering Lab 3 <input type="checkbox"/> Cell Engineering Lab 4 (※)		
<input type="checkbox"/> Pathogenic viruses (<input type="checkbox"/> New · <input type="checkbox"/> Renewal)	Approved number :		Supervisor's name : (Affiliation :)	
	Approved room	<input type="checkbox"/> Pathogenic Virus Lab (※) <input type="checkbox"/> Cell Engineering Lab 4 (※)		
<input type="checkbox"/> Animal experiments (<input type="checkbox"/> New · <input type="checkbox"/> Renewal)	Approved number :		Supervisor's name : (Affiliation :)	
	Approved room	<input type="checkbox"/> X-ray Fluoroscopy Lab (※) <input type="checkbox"/> X-ray Irradiation Lab <input type="checkbox"/> Wet Lab <input type="checkbox"/> 3T MRI Lab <input type="checkbox"/> Animal Imaging Lab <input type="checkbox"/> Laser Microscope Lab 1 <input type="checkbox"/> Cell Engineering Lab 3		

【NOTE】 Any person who is engaged in the operation of recombinant DNA experiments, experiments handling pathogenic microorganisms, animal experiments must comply with the regulations in life science.