Central Research Laboratory

2025 Guidance for New Users Research Equipment Section Central Research Laboratory

Central Research Laboratory
Shiga University of Medical Science

2025 Guidance for New Users

Research Equipment Section, Central Research Laboratory (CRL)

1)	Overview: Usage of Research Equipment Section	1
2)	Industrial safety and health	2
3)	How to use the 'recombinant DNA room' and 'pathogenic virus room'	4
4)	How to start animal experiments at CRL	6
5)	Research support by CRL	7
6)	How to use X-ray generators at CRL	8
7)	Introduction of each equipment available at CRL	9
8)	How to avoid unnecessary accidents?	13

Appendix: Equipment List / CRL Floor Map / Staff List

***** How to Register *****

* After taking the guidance, fill out the "Application form for CRL User Registration" referring to the "Example of Entry". Submission of the form will complete your registration.

* Once the registration has been processed;

* Registered users can continue their registration by simply going through the renewal procedure (submitting the renewal application form) from the following year.

Central Research Laboratory, SUMS
Tel: 077-548-2300
E-mail: hqcrl@belle.shiga-med.ac.jp
URL: <http://wwwcrl.shiga-med.ac.jp/index.html>

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Guidance for New Users FY2025 Research Equipment Section Central Research Laboratory

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1) Overview: Usage of Research Equipment Section

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For on site attendee

1. Attend this lecture
2. Download "Application Form" and submit it by e-mail attachment
(hqcri@belle.shiga-med.ac.jp)

For Zoom on line attendee

1. Attend this lecture via Zoom
2. We will provide a unique link in the Zoom chat during the lecture
3. Go to the link and download both "Survey form" and "Application Form"
4. Fill out both forms and send them by e-mail attachment for approval

* All users need to submit a renew application every year

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1. University level

- ・All members need to comply the Industrial Safety and Health Act
- ・All members need to take the education program regarding safety management

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1. Make your reservation on the notebook placed near the equipment; Write your name, affiliation, and telephone number on your desired date/time slot.
2. Do not make reservation long time. The reservation will be canceled if no one use it 30 minutes after the start time.
3. Need to bring consumables necessary for using the equipment.
4. Follow the designated rule for each equipment.
5. Fill out the notebook after use (This is important for calculating the user charge).
6. Clean up the equipment.

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Follow the designated rule for each equipment

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- ・When you use the equipment first time, **you must ask the CRL staffs and learn** how to use it.

- ・User manuals are available on the CRL website
http://www.crl.shiga-med.ac.jp/home/kiki_bumon/g_book/contents.html

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-
1. The CRL facility is open 24 hours.
Asks the CRL staff if you need to use locked rooms.
 2. No eating, drinking, and smoking allowed.
 3. Additional registration is necessary for access to 'Genetic Engineering room', 'P2 Pathogenic Virus room', and 'Animal X-ray room'.
 4. CRL provides "Special service" for several instruments and analyzers.

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Participants in the above seminars can obtain credits for the PhD course lectures.

Details are indicated at the CRL website.
(<http://www.crl.shiga-med.ac.jp/home/seminar/seminar.html>)

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Formaldehyde should be used under the fume hood in the CRL.
Users should sign in a Formaldehyde Log after use.

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Evacuate from the fire zone if you think it dangerous.

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
How to use the "recombinant DNA room" and "pathogenic virus room"

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Act on the Conservation and Sustainable Use of Biological Diversity through Regulations on the Use of Living Modified Organisms(LMOs), so-called **Cartagena Law**

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Containment measures (Type 2 Use of LMOs)
Possible uses are those with containment measures to prevent to dispersal of LMOs in the environment (Uses in laboratories etc).

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
Containment Measures of P1 Level Laboratory

Closed doors and windows on P1 experiment
Limited access
Washed hands
Inactivated all contaminated things

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Containment Measures of P1AL (experiment with *Animal*) Level Laboratory

遺伝子組換え生物等が逃亡その他拡散しない構造の容器に入れること
When a LMO is taken out of the laboratory in process of an experiment, the LMO shall be put in a container of the structure that prevents it from leaking or other dispersion.

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Keep the recombinant DNA room clean

初めて使用する場合、必ずセンター職員に聞く
When using the equipment for the first time, please ask the CRL staffs how to use it.

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Containment Measures of P2 Level Laboratory

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Application Process for Use of CRL Pathogenic Virus room

```

graph TD
    A[支援センター機器部門利用登録  
User Registration in Research Equipment Section, CRL  
http://gakuorl.shiga-med.ac.jp/gakunai/gakunai/0_riyousinsei/kiki_bumon/sinki.html] --> B[滋賀医科大学バイオセーフティ委員会への届出  
Application to the SUMS Biosafety Committee  
http://gakuorl.shiga-med.ac.jp/gakunai/gakunai/0_riyousinsei/kiki_bumon/virus2.html]
    B --> C[実験実習支援センター病原体等取扱申込書の提出  
Submission of the "CRL Application Form for Handling Pathogens, etc."  
http://gakuorl.shiga-med.ac.jp/gakunai/gakunai/6_Regulations/kiki_regulations/kiki_regulations_virus_0704.html]
    C --> D[支援センター管理責任者の承認  
Approval by the CRL Director]
    D --> E[病原ウイルス実験室の利用  
Acceptance for Use of Pathogenic Virus Lab]
    
```

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4) How to start animal experiments at CRL

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Need to complete following:
1. Completed the training program in animal experiment.
2. Pass the qualification exam of animal experiment.
3. Obtain the certification on animal experiment.

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Animal experiments in the CRL rooms

「動物の愛護及び管理に関する法律」、「滋賀医科大学動物実験規程」に従う。
Perform animal experiments under the Law Concerning the Protection and Control of Animals and the regulations of Research Center for Animal Life Science in SUMS.

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
You have 2 choices

- 1) prepare sequencing reaction and submit the sample for running
- 2) submit a plasmid for DNA sequencing.

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
自動包埋装置
Automate tissue Processor

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Yamamoto (Ext. 2304)
4F Room401

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Additional support provided by CRL

Bioinformatics解析（門田先生）

Bioinformatics analysis by Dr. Kadota

10x Genomics Single cell RNA-sequencing reagentの代理購入

Provide 10x Genomics Single cell RNA-sequencing reagent

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Provide qPCR plates

MOLSIS Molecular Simulation & Informatics Systems Software

FlowJo software for Flow cytometer and FACS analysis

Leica LAS-X software for image processing

Photoshop

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Hands-on guidance on your experiments

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A young researcher in the SUMS hospital wanted to determine gene mutations in the patient sample. The CRL staff helped designing primers, PCR, and DNA sequencing. The mutation was successfully detected.

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Please ask us if you would like to know how to perform experiments. We may not help your issues directly, but will help your research as much as we can.

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Technical staff, Central Research Laboratory

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Ordinance and regulations of X-ray generators

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Ordinance on Prevention of Ionizing Radiation Hazards

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Basic Principle: Minimize the exposure of workers to ionizing radiation, as low as possible

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SUMS Regulations on the Regulation of Radioisotopes, etc.

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CRL Rules on X-ray Generator Utilization

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X-ray Irradiator
(1st floor, Research Center for Animal Life Science)

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X-ray Diffractometer
(1st floor, Central Research Laboratory)

No risk of x-ray leaking from these generators

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Registration of use of X-ray generators

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Take an education and training class (1 hour) provided by CRL

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Registration is completed for users of X-ray irradiators
(No need blood test and wearing a film badge)

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Users need to wear a pocket dosimeter

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X-ray user must record the radiation dose to the logbook.

X-ray fluoroscope for animals

Users should use appropriate protective wears to prevent inevitable exposure of X-ray.

Laboratory Animal X-ray Fluoroscope in Central Research Laboratory

Users will be exposed to X-ray!

Registration of use of X-ray fluoroscope for animals

Registration

●血液検査
To use X-ray fluoroscope, users need to take blood tests

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Registration as a radiation worker in SUMS

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担当: 山元(内線2304)
Technical staff: T. Yamamoto (Ext. 2304),

Introduction of each equipment available at CRL

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For solid samples

For liquid samples

NMR can determine the molecular structure of a sample placed in a strong magnetic field by nuclear magnetic resonance.

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The device is capable of measuring molecular weight of small molecules.

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This device can measure molecular weight of gas molecules.

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HPLC can separate and quantify components dissolved in a liquid sample.

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Peptides and proteins can be quantified and purified.

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
Compatible to 6-384 well Plates for measuring absorption, fluorescence, and luminescence.

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Nucleic acid, protein, etc. can be quantified with a small sample of 1 - 2 μ L.

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Quantify nucleic acids or proteins using a microchip

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
For quantitative PCR analysis of gene expression. Users may purchase PCR plates from CRL.

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Chemiluminescence detection for Western blotting.

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Detection of nucleic acid stained with EtBr by LED

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Ask the CRL staff which rotor will be compatible.

Ultracentrifuge

CRL 3rd floor (To be relocated to RCALS)

In vivo imaging system NEWTON 7.0

Detection of bioluminescence and fluorescence in live animals

New

High-resolution Ultrasonic Echo Imaging System

Visualize *in vivo* anatomical, functional, and physiological data in real time using ultrasound echo.

CRL 3rd floor

Redox Analyzer REDOXLIBRA

Measurement ROS and free radicals

Chemistry Analyzer VetScan

Measurement blood chemistry

Blood Analyzer Celltac α

Measurement of blood cell parameters

New

Tetrapolar Bioimpedance Spectroscopy Device ImpediVET (To be relocated to RCALS)

Measurement of fluid status and tissue composition in live animals

CRL 3rd floor

Ultrasonic cell crusher

Homogenization with beads Minilys

Polytron homogenizer

Vacuum freeze dryer

CRL 4th floor

Confocal Laser Scanning Microscope

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
Confocal laser scanning microscopy is a technique for obtaining high-resolution optical images with depth selectivity.

Please ask CRL staff which microscope will be suitable for your experiment.

CRL 4th floor

Other microscopes

1. Inverted Fluorescence Microscope IX83 (Olympus)
2. DeltaVision High-resolution Fluorescence Microscope (GE)
3. All-in-one fluorescence microscope BZ9000 (Keyence)
4. Eclipse Ni-E (Nikon)

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Ask us which microscope will be suitable for your experiment.

CRL 4th floor

New


Transmission Electron Microscope HT7800 (Hitachi)

Scanning Electron Microscope JSM-7505FA (JEOL) JSM-6010LA (JEOL)

CRL 4th floor

Ultramicrotome ULTRACUT E (Reichert-Jung)
Making thin sections for electron microscope

Draft Chamber
Use organic solvent under the draft chamber


CRL 4th floor

Automated Tissue Processor
Paraffin embedding

Tissue Embedding Console
Making paraffin blocks

Microtome
Making paraffin sections


Cryostat
Making frozen sections

CRL 4th floor

Flow Cytometer

1. LSRFortessa X-20 (BD)
2. FACSCanto II (BD)
3. CytoFLEX S (Beckman)
4. FACS Calibur (BD)
5. Attune (AB)

New

Cell analysis

FACS

1. FACS Aria Fusion (BD)
2. FACS Aria Fusion (BD)

Cell analysis and sorting

Ask us which equipment will be suitable for your experiment.

CRL 4th floor

autoMACS Separator (Miltenyi Biotec)

Cell separation using antibodies

Suspension Array System Bio-Plex200 (Bio-Rad)

Measuring cytokines using special beads

Extracellular Flux Analyzer XFe24 (Seahorse Bioscience)

Measuring metabolism in cultured cells

CRL 4th floor

CO₂ Incubator and clean bench

Users may request use of CO₂ incubator or clean bench for cell culture.

CRL 4th floor

Large Capacity Liquid Nitrogen Cryopreservation Equipment (MVE Heco1536P)
Store cells in liquid nitrogen gas

Liquid Nitrogen
Users can buy liquid nitrogen from CRL

Milli-Q and Elix Water System (2nd and 4th floor)
Users can buy pure water from CRL

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8) How to avoid unnecessary accidents?

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The rotor came off from the spindle and rolled around in the chamber. The rotor and spindle were severely damaged.

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The bucket came off from the running rotor at 40,000 rpm and crashed onto the chamber wall at a velocity of over 1,600 km/h. Relative Centrifugal Field is 193,000 g.

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Tips for safety centrifuge operation

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
Leica TCS SP8 X
Confocal Laser Scanning Microscope

In 2023, someone broke the objective lens. It cost 700,000 yen for replacement.

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Nikon/AX Confocal Microscope

The user wanted to see the edge of the slide and moved the objective lens upward to adjust focus. Then, the objective lens hit and bent the stage holder.

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-
1. If you cause any accident, immediately inform it to the CRL staff.
 2. If you notice anything unusual such as a trouble with equipment, immediately inform it to the CRL staff.
 3. If you are not sure how to use the equipment, ask the CRL staff about it.
DO NOT use the equipment only on your own idea and judgment.
 4. In using the CRL facilities, the rules at CRL are prior to those at your lab.
 5. DO NOT follow the user who acts in incorrect manner and method.

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On site attendee

After this lecture, download the application form and submit it to CRL via e-mail

On line Zoom attendee

We will display a link for download site of two forms to be submitted ("survey form" and "application form for CRL User Registration") in the Zoom chat box during this course. By sending them using email, your registration will be completed.

Contact Number & E-mail Address

Research Equipment Section of Central Research Laboratory

○ Director of CRL

Job title	Name	Ext.	E-mail
Professor	Itoh, Yasushi	2171	yasushii@belle.shiga-med.ac.jp

○ Faculty of CRL

Job title	Name	Ext.	E-mail
Associate Professor	ASAHINA, Kinji	2301	asahina@belle.shiga-med.ac.jp
Senior Associate Professor	TOYODA, Futoshi	2322	toyoda@belle.shiga-med.ac.jp

○ Technical Staff

Job title	Name	Ext.	E-mail
Technical Staff	MORI, Yasuhiro	2307	moriyasu@belle.shiga-med.ac.jp
Technical Staff	OKAMOTO, Kumi	2302	okamotok@belle.shiga-med.ac.jp
Technical Staff	FUKUNAGA, Sachiko	2303	sfuku@belle.shiga-med.ac.jp
Technical Staff	TERADO, Tokio	2306	terado@belle.shiga-med.ac.jp
Technical Staff	YAMAMOTO, Takefumi	2304	tyama@belle.shiga-med.ac.jp

○ Administrative Staff

Job title	Name	Ext.	E-mail
Technical Staff	YAMAKAWA, Nobuko	2300	nyama@belle.shiga-med.ac.jp

○ Central Research Laboratory

Location	Ext.	E-mail
2F Central Research Laboratory office	2300	hqcrl@belle.shiga-med.ac.jp

More information on CRL is available at the CRL web site (<http://www.crl.shiga-med.ac.jp/index.html>) where you can get a list of installed equipment, online edition of the operation manuals, etc. and download various request/application forms.

Device List (Research Equipment Section, Central Research Laboratory)

BIOCHEMISTRY	DEVICE/SYSTEM	FEE (JPY)	LOCATION
Nuclear Magnetic Resonance Spectroscopy	JEOL/JNM-ECZ400S, JNM-ECX400	1,000/h (Day Max 5,000)	1F NMR Lab
Microplate Reader	TECAN/Infinite M PLEX	100/use	2F Biochemical Analysis Lab 1
Spectrophotometer	Thermo Fisher/NanoDrop One, NanoDrop ND-1000	—	2F Molecular Biology Analysis Lab
UV-Vis-NIR Spectrophotometer	JASCO/V-570DS	—	2F Biochemical Analysis Lab 1
Fourier Transform Infrared Spectrometer	SHIMADZU/IRPrestige-21	—	2F Protein Engineering Lab
Fluorescence Spectrometer	HITACHI/F-2500	—	2F Protein Engineering Lab
Circular Dichroism Spectrometer	JASCO/J-1500	—	2F Biochemical Analysis Lab 1
Gel Imaging System	NIPPON Genetics/FAS5 UVP/GDS-7900	—	2F Molecular Biology Analysis Lab
Lumino Image Analyzer	FUJIFILM/LAS-4000	100/use	2F Molecular Biology Analysis Lab
Chemiluminescence Imaging System	Vilber-Lourmat/FUSION	100/use	2F Molecular Biology Analysis Lab
Preparative Ultracentrifuge	Beckman/Optima L-90K	100/10min (5,000/More than 500min)	2F Centrifuge Room
Tabletop Ultracentrifuge	Beckman/Optima MAX-TL	100/10min (5,000/More than 500min)	2F Centrifuge Room
High Speed Refrigerated Centrifuge	Beckman/HP-25·HP-26	200/use	2F Centrifuge Room
Hybrid High Speed Refrigerated Centrifuge	KUBOTA/6200	—	2F Centrifuge Room
Low Speed Centrifuge	TOMY/LC-131	—	2F Molecular Biology Analysis Lab
High Speed Refrigerated Centrifuge (P2 Lab)	Beckman/Avanti J-E	200/use	2F Gene Engineering Lab P2
High Speed Refrigerated Microcentrifuge	TOMY/MX-305	—	2F Gene Engineering Lab P2
	HITACHI/CF16RN	—	3F Centrifuge Room
Plate Centrifuge	KUBOTA/PlateSpin II	—	2F Molecular Biology Analysis Lab
Multipurpose Centrifuge	Beckman/GS-15R	—	2F Gene Engineering Lab P2
Low Speed Refrigerated Centrifuge	Himac/CF7D2	—	2F MS Spectrometry Lab
Vacume Concentrator/Speedvac Concentrator	Thermo Fisher/SPD1010	200/use	2F Centrifuge Room
Shaking Incubator	Sanki Seiki/SCS-R·SCS-12R, TAITEC/BR-30L, BR-30LF, BR-43FL	—	2F Gene Engineering Lab P2
Lyophilizer	LABCONCO/LL-1	200/d	3F Cell Homogenization Room
Ultrasonic Homogenizer	BRANSON/Sonifier Model 450- Advanced	—	3F Cell Homogenization Room
Cap Horn Homogenizer	MISONIX/ASTRASON MODEL XL2020	—	3F Cell Homogenization Room
Polytron Homogenizer	KINEMATICA/MODEL PT1200E	—	3F Cell Homogenization Room
Bead Mill Homogenizer	Bertin Technologies/Minilys	100/use	3F Biological Sample Analysis Lab 1

BIOCHEMISTRY	DEVICE/SYSTEM	FEE (JPY)	LOCATION
Chromatography Chamber	NIHON FREEZER/MC-30EF3	500/rack/m	3F Frozen Strage Room
Electronic Balance	SHIMADZU/AUW320	—	2F Molecular Biology Analysis Lab
High Performance Liquid Chromatograph	SHIMADZU/Nexera	unfixed	2F MS Spectrometry Lab
	SHIMADZU/HPLC System	200/use	2F Biochemical Analysis Lab 1
HPLC AKTA System	GE Healthcare/AKTA purifier 100	200/use	2F Biochemical Analysis Lab 1
Gas Chromatograph Mass Spectrometer	SHIMADZU/GCMS-QP5050A	100/use	2F MS Spectrometry Lab
Time-of-Flight Mass Spectrometer	SHIMADZU/MALDI-TOF/MS AXIMA Confidence	1,000/use	2F MS Spectrometry Lab
X-ray Diffractometer	Philips/MPD-1880	1,000/use	1F X-ray Irradiation Lab
DNA Sequencer	ABI/PRISM 3130XL GeneticAnalyzer	(ELP)200/sample (ELP+Seq)600/sample	3F Gene Analysis Lab
	HITACHI/DS3000	(ELP)250/sample (ELP+Seq)650/sample	3F Gene Analysis Lab
Acoustic Solubilizer	Covaris/S220	100/use	2F Molecular Biology Analysis Lab
Real Time Quantitative PCR (LightCycler)	Roche/LightCycler480 System II	100/use (Plate)8,052/set	2F Molecular Biology Analysis Lab 2F CRL Office
Thermal Cycler	ABI/GeneAmp PCR System 2720, Roche/LightCycler, Thermo /MiniAmp Plus	—	2F Molecular Biology Analysis Lab
Electrophoresis Apparatus	Pharmacia/IPG Phor	500/use	2F Protein Engineering Lab
Bioanalyzer	Agilent Technologies/Agilent 2100 Bioanalyzer Limited	100/use	2F Molecular Biology Analysis Lab
Suspension Array System	Bio-Rad/Bio-Plex200	3,000/use	4F Cell Engineering Lab 4
Sample Derivatization System	Thermo Fisher/Reacti-Therm	200/h	2F Protein Engineering Lab
Sample Concentrator		200/h	2F Protein Engineering Lab
MOE (Molecular Operating Environment)	MOLISIS	—	3F Gene Analysis Lab

MORPHOLOGY	DEVICE/SYSTEM	FEE (JPY)	LOCATION
Transmission Electron Microscope	HITACHI/HT7800	unfixed	4F Electron Microscope Lab
Field Emission Scanning Electron Microscope	JEOL/JSM-7505FA	2,000/h	4F Electron Microscope Lab
Scanning Electron Microscope	JEOL/JSM-6010LA	2,000/h	4F Electron Microscope Lab
Ultramicrotome	Reichert/ULTRACUTE	500/use	4F Microscopy Sample Preparation Room 1
Color CCD Camera	Nikon/ECLIPSE Ni-E	100/use	4F Image Processing Room
	NIPPON ROPER/Retiga2000R	100/use	4F Cell Culture Room
Live Cell Imaging Microscope	Andor/Dragonfly201	300/h(Day Max 5,000)	4F Confocal Laser Scanning Microscope Lab 1
High-resolution Fluorescence Microscope	GE/DeltaVision Elite	200/use	4F Microscope Lab 1
All-in-One Fluorescence Microscope	KEYENCE/BIOREVO BZ-9000	2,000/h	4F Image Processing Room

MORPHOLOGY	DEVICE/SYSTEM	FEE (JPY)	LOCATION
Inverted Fluorescence Microscope	OLYMPUS/IX83	100/use	4F Microscope Lab 2
Confocal Laser Scanning Microscope	LEICA/STELLARIS 8	1,000/use	4F Confocal Laser Scanning Microscope Lab 2
	LEICA/TCS SP8 X	1,000/use	4F Confocal Laser Scanning Microscope Lab 3
	Nikon/AX with NSPARC	1,000/use	4F Fluorescence Microscope Lab
	OLYMPUS/FV1000-D	1,000/use	3F Biological Sample Analysis Lab 1
Microtome	LEICA/SM2010R	100/use	4F Microscopy Sample Preparation Room 2
Tissue Processor	SAKURA/Tissue-Tek VIP 6-J0, etc.	1,000/use	4F Microscopy Sample Preparation Room 2
Tissue Embedding Console	SAKURA/Tissue-Tek 6	100/piece	4F Microscopy Sample Preparation Room 2
Cryostat	LEICA/CM3050 S	500/use	4F Microscopy Sample Preparation Room 2
3D/4D Image Analysis and Visualization Software	BITPLANE/Imaris	100/use	4F Image Processing Room

CELL CULTIVATION	DEVICE/SYSTEM	FEE (JPY)	LOCATION
Nucleofector	LONZA/4D-Nucleofector, amaxa/Nucleofector	100/use	4F Cell Culture Lab
CO2 Incubator	YAMATO/IP400, SANYO/MCD-175, ESPEC/BNA-111, HIRASAWA/CPD-2701	45,000/y + CO2 cost	4F Cell Culture Lab
Clean Bench (Laminar Flow Cabinet)	ESCO/Airstream PRO ESC-BCB-4A7, Showa Kagaku/S-1801WBV, etc.	—	4F Cell Culture Lab
Flexercell	FLEXERCELL/FLEX-I	2,000/use	4F Cell Culture Lab
Microinjection System	Eppendorf/InjectMan NI2, OLYMPUS/IX-73	100/use	4F Cell Engineering Lab 3
Flow Cytometer	BD/FACSCalibur·FACSCanto II	1,500/use	4F Cell Engineering Lab 2
	BD/LSRFortessa	1,500/use	4F Cell Engineering Lab 4
	Beckman Coulter/CytoFLEX S	1,500/use	4F Cell Engineering Lab 1
Cell Sorter	BD/FACSAria Fusion	5,000/use (Analysis)1,500/use	4F Cell Engineering Lab 4
Magnetic Cell Sorter	Miltenyi Biotec/autoMACS Pro Separator	3,500/use (No consumables)1,100/use	4F Cell Engineering Lab 1
Cell Dynamic Analysis System	GE/EZ-TAXIScan	unfixed	4F Cell Culture Lab
Extracellular Flux Analyzer	Seahorse/XFe24	100/use	4F Cell Culture Lab
Deep Freezer (Frozen Strage Room) Container Basket	REVCO/ULT-1490-3J-D30, etc.	(Container)6,000/y (Drawer Case)4,000/y	3F Frozen Strage Room
Large Liquid Nitrogen Sample Storage Container	CHART/MVE Heco1536P-190AF·GB	550/m	4F Workroom
Biosafety Cabinet	ESCO/ClassII Type A2	—	4F Cell Engineering Lab 4
	Thermo/1387	—	4F Pathogenic Virus Lab
Preparative Ultracentrifuge	Beckman/Optima L-90K	100/10 min (5,000/More than 500min)	4F Pathogenic Virus Lab
High Speed Refrigerated Centrifuge	Beckman/Avanti J-E	200/use	4F Pathogenic Virus Lab


CELL CULTIVATION	DEVICE/SYSTEM	FEE (JPY)	LOCATION
Tabletop Centrifuge	KUBOTA/KS-5000	—	4F Pathogenic Virus Lab
Tabletop Refrigerated Centrifuge	Beckman/Allegrax-30R	—	4F Cell Culture Lab
Ethylene Oxide Gas Sterilizer	ELK CORPORATION/SA-360ECO IKI/FRH36	3,000/use	4F Frozen Storage Room
X-ray Irradiator	AcroBio/RS320	1,000/use	Research Center for Animal Life Science

PHYSIOLOGY / ANIMAL	DEVICE/SYSTEM	FEE (JPY)	LOCATION
Microdialysis	EICOM/HTEC-500	1,000/use	2F Biochemical Analysis Lab 1
Animal Blood Gas Analyzer	ABAXIS/VetScan VS2	—	3F Biological Sample Analysis Lab 1
Automatic Blood Cell Counter	NIHON KOHDEN/ Celltac α MEK-6550	200/sample	3F Biological Sample Analysis Lab 1
Oxidation Redox Analyzer	WISMERLL/REDOXLIBRA	100/use	3F Biological Sample Analysis Lab 1
Tetrapolar Bioimpedance Spectroscopy (BIS) Device	Bio Research Center/ImpediVET	100/use	3F Biological Sample Analysis Lab 1
C-arm X-ray Fluoroscope for Laboratory Animals	SIEMENS/ARCADIS Avantic Gen2	1,000/use	1F Animal X-ray Fluoroscopy Lab
Heart Perfusion System	ADInstruments/Langendorff	1,000/use	3F Biological Sample Analysis Lab 1
High-resolution Ultrasonic Echo Imaging System for Small Animals	VisualSonics/Vevo2100	250/15min (Day Max 5,000)	3F Small Animal Imaging Analysis Lab
In vivo Imaging System for Small Animals	Vilber-Lourmat/NEWTON 7.0	200/h	3F Small Animal Imaging Analysis Lab
3-Tesla MRI Scanner	SIEMENS Healthineers/ MAGNETOM Verio dot 3T	400,000/y 5,000/d, 1,500/h	1F 3T MRI Lab
4.7-Tesla MRI Scanner	BRUKER/BioSpec47/40USR	400,000/y 5,000/d, 1,500/h	Research Center for Animal Life Science


ADDITIONAL SERVICE	DEVICE/SYSTEM	FEE (JPY)	LOCATION
Elix Water	MILLIPORE/Milli-Q Integral 5	10/L	2F Biochemical Analysis Lab 4F Workroom
Milli-Q Water	MILLIPORE/Milli-Q Integral 5	80/L	2F Biochemical Analysis Lab 4F Workroom
Ice Machine	Scotsman/F0522B Hoshizaki/FM-340AK-SA	—	2F Biochemical Analysis Lab 4F Workroom
Liquid Nitrogen	100 liter, 50 liter tank	450/L (Market Price)	4F Workroom
Oxygen Gas	purity 99.5%, 4.0Kg/cm ²	4,246/cylinder (Market Price)	1F Bombe Room 1F Wet Lab
Compressed Air		5,500/cylinder (Market Price)	1F Wet Lab
Human Sample Storing	REVCO, NIHON FREEZER	2,000/freezer/m	1F Human Sample Resource Room

WET LAB	DEVICE/SYSTEM	FEE (JPY)	LOCATION
Electric Knife	ellman/SURGITRON Dual EMC 90	—	1F Wet Lab
Small Animal Anesthetizer	SHIN-EI INDUSTRIES/A.D.S.1000	—	1F Wet Lab

1F·2F of Central Research Laboratory

3F·4F of Central Research Laboratory


Application Procedures for New User Registration

Those who wish to use the facilities and equipment of Central Research Laboratory (CRL) must be registered users. To be the user, they are required to take the Guidance for New Users managed by CRL. Depending on the nature of the experiments they plan to conduct, other application procedures may be required in addition to the user registration.

They must complete the necessary procedures according to the section they wish to use and the type of experiments they are planning to conduct.

Once you have completed the registration ...

Each piece of equipment at the CRL has its own rules on how to reserve and use it. If you use equipment for the first time, please be sure to ask the staff member in charge of the equipment how to use it before using it, rather than having your lab members teach you how to use it.

Items that Require Application to the University

- ✓ Recombinant DNA Experiments: Application to the SUMS Genetic Recombination Experimentation Safety Committee (Research Promotion Division)
- ✓ Experiments Handling pathogenic microorganisms: Application to the SUMS Biosafety Committee (Research Promotion Division)
- ✓ Experimental Animals: Application to the SUMS Animal Experimentation Committee (Research Center for Animal Life Science)
- ✓ Experiments with Animal X-ray Fluoroscope: Registration as an X-ray photographing worker (Research Promotion Division)
- ✓ Experiments with X-ray Irradiator: Registration as an X-ray worker (Research Promotion Division)
- ✓ Human Subjects Research: Application to the SUMS Ethics Committee

Instructions and Notes for Completing the “Application Form for CRL User Registration”

Those who conduct the following experiments at CRL must enter the approval number and other necessary information on the Form, and attach the prescribed documents, if necessary:

- (1) **Recombinant DNA experiments:** Those who performing P1 and/or P2 level experiments must include their approval number. Additionally, the Gene Engineering Lab and/or Cell Engineering Lab users must contact the staff (Ext.2306) to register their ID card information for the entry/exit managing system.
- (2) **Experiments handling pathogenic microorganisms:** Those who performing the experiments of this kind must include their approval number and submit the prescribed documents. Additionally, the Pathogenic Virus Lab and/or Cell Engineering Lab #4 users must contact the staff (Ext.2307) to register their ID card information for the entry/exit managing system.
- (3) **Animal experiments:** Those who performing the animal experiments must include their approval number.
- (4) **X-ray experiments** (Animal X-ray Fluoroscope, X-ray Irradiator, X-ray Diffractometer): X-ray generator users are required to take the X-ray training course prior to using the equipment.
- (5) **Human Sample Resource Room:** The HSR users must submit the prescribed form. Additionally, they must register their ID card information for the entry/exit managing system.

【Caution】 Matters to be entered, documents to be submitted

Items	Approval number	Entry/exit managing system	Document to submit	Documents to submit	Staff (Ext.)
Recombinant DNA experiments	●	●			Terado (2306)
Handling pathogenic microorganisms	●	●	●	①CRL Application Form for Handling Pathogens, ②Application for Handling Pathogens (Form 1 or 2)	Mori (2307)
Experimental Animals	●				Terado (2306)
Animal X-ray Fluoroscope		●			Yamamoto (2304)
Wet Lab	●				Terado (2306)
Human Sample Resource Room		●	●	HSR Access Application	Yamamoto (2304)

EXAMPLE

Take extra care not to make mistake or omit anything

2025 New User Registration Form

Central Research Laboratory (CRL), Research Equipment Section

I will abide by rules and regulations established by CRL.

Date	2025 / 00 / 00	User No. (CRL Entry)	
Affiliation	Department / Center of 00		
Position	Graduate student		
Name	SHIGA Taro		
Staff (Student) ID number	88888888		
Extension/PHS number	2300		
e-mail address	hqcl	@belle.shiga-med.ac.jp @g.shiga-med.ac.jp	
Head of department			

Your research projects	<input checked="" type="checkbox"/> Check the equipment that you plan to use		
Items you plan to use	<input type="checkbox"/> Spectrophotometer <input type="checkbox"/> Mass spectrometer <input type="checkbox"/> Nuclear magnetic resonance <input type="checkbox"/> HPLC <input type="checkbox"/> Temperature control centrifuge <input type="checkbox"/> Ultracentrifuge <input type="checkbox"/> Lyophilizer <input type="checkbox"/> Vacuum concentrator <input type="checkbox"/> Micro-dialysis analyzer <input type="checkbox"/> Homogenizer <input type="checkbox"/> DNA sequencer <input type="checkbox"/> Bioanalyzer <input type="checkbox"/> Covaris solubilizer <input type="checkbox"/> Bio-Plex suspension array <input type="checkbox"/> MOE (Molecular Operating Environment) <input type="checkbox"/> Chemilumi imaging <input type="checkbox"/> PCR/qPCR	<input type="checkbox"/> Gel imaging <input type="checkbox"/> Chromatography chamber <input type="checkbox"/> Electron microscope <input type="checkbox"/> EM sample preparation <input type="checkbox"/> Light microscope <input type="checkbox"/> Microscope sample preparation <input type="checkbox"/> Image analyzer <input type="checkbox"/> Laser microscope <input type="checkbox"/> Live imaging microscope <input type="checkbox"/> Fluorescence microscope <input type="checkbox"/> CO ₂ incubator <input type="checkbox"/> MACS <input type="checkbox"/> Flow cytometer <input type="checkbox"/> Cell sorter <input type="checkbox"/> EZ-TAXIScan <input type="checkbox"/> Nucleofector <input type="checkbox"/> Flux analyzer	<input type="checkbox"/> Deep freezer <input type="checkbox"/> Liquid nitrogen dewar <input type="checkbox"/> Autoclave <input type="checkbox"/> Heart perfusion system <input type="checkbox"/> Small animal imaging <input type="checkbox"/> in vivo imaging <input type="checkbox"/> VetScan <input type="checkbox"/> X ray fluoroscope <input type="checkbox"/> X ray generator <input type="checkbox"/> MRI scanner <input type="checkbox"/> Millipore water <input type="checkbox"/> Liquid nitrogen <input type="checkbox"/> Oxygen gas <input type="checkbox"/> Handicraft tools <input type="checkbox"/> Human sample storage <input type="checkbox"/> Wet lab

Please fill in the following fields only if your experiment at CRL is involved in recombinant DNA, pathogenic viruses and microorganisms, or animal use.

※ Access control laboratories with ID card authentication

<input checked="" type="checkbox"/> Recombinant DNA <input type="checkbox"/> New <input type="checkbox"/> Renewal	Approved No.	5-1	Supervisor's name		Affiliation	
	Approved room	<input checked="" type="checkbox"/> Genetic Engineering Lab ※ <input type="checkbox"/> 3T MRI Lab <input type="checkbox"/> Biological Sample Analysis Lab 1 <input type="checkbox"/> Cell Engineering Lab 1 <input type="checkbox"/> Cell Engineering Lab 3				
<input checked="" type="checkbox"/> Pathogenic viruses <input type="checkbox"/> New <input type="checkbox"/> Renewal	Approved No.	5-1	Supervisor's name		Affiliation	
	Approved room	<input checked="" type="checkbox"/> Pathogenic Virus Lab ※ <input type="checkbox"/> Cell Engineering Lab				
<input checked="" type="checkbox"/> Animal experiments <input type="checkbox"/> New <input type="checkbox"/> Renewal	Approved No.	2024-1	Supervisor's name		Affiliation	
	Approved room	<input checked="" type="checkbox"/> X-ray Fluoroscopy Lab ※ <input type="checkbox"/> 3T MRI Lab <input type="checkbox"/> Biological Sample Analysis Lab 1				

- ☐ X-ray Irradiation Lab
☐ Animal Imaging Lab
☐ Cell Culture Lab
☐ Cell Engineering Lab 2
☐ Cell Engineering Lab 4

contact the staff (ext. 2307) and follow the prescribed

- ☐ Wet Lab
☐ Animal Imaging Lab
☐ Cell Engineering Lab 3

[NOTE] Any person who is engaged in the operation of recombinant DNA experiments, experiments handling pathogenic microorganisms, animal experiments must comply with the regulations in life science, research.