Central Research Laboratory

2023 Guidance for New Users Research Equipment Section Central Research Laboratory

Central Research Laboratory
Shiga University of Medical Science

2023 Guidance for New Users
Research Equipment Section, Central Research Laboratory (CRL)

1)	Overview: Usage of Research Equipment Section, CRL	1
2)	Industrial safety and health	4
3)	How to use the recombinant DNA laboratory	7
4)	How to use experimental animals at CRL	10
5)	Special service rendered by CRL	12
6)	Procedures for using X-ray generators	14
7)	Equipment introduction #1 (Instruments on the 1st and 4th floors)	16
8)	Equipment introduction #2 (Instruments on the 2nd and 3rd floors)	21

Appendix: Equipment List / CRL Floor Map / Staff List

```
* ***** How to Register ***** *
*
* After taking the guidance, you can download the "Application form for CRL
* User Registration". By sending it using email, your registration will be
* completed.
*
* Once the registration has been processed...
* Registered users can continue their registration by simply going through
* the renewal procedure (submitting the renewal application form) from the
* following year.
* *****
```

Central Research Laboratory, SUMS
Tel: 077-548-2300
E-mail: hqcrl@belle.shiga-med.ac.jp
URL: <http://wwwcrl.shiga-med.ac.jp/index.html>

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Guidance for New Users FY2023 Research Equipment Section, Central Research Laboratory

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1) Overview: Usage of Research Equipment Section, CRL

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User registration process at Research Equipment Section, CRL

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For on site attendee

1. Attend this lecture
2. Download "Application Form" and submit it by e-mail attachment
(hqrcr@belle.shiga-med.ac.jp)

For Zoom on line attendee

1. Attend this online Zoom guideline lecture
2. We will provide a unique link in the Zoom chat during the lecture
3. Go to the link and download both "Survey form" and "Application Form"
4. Fill out both forms and send them by e-mail attachment for approval

* By annual renew application, the user registration is maintained.

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Why is the user registration required?

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1. Request from National University Corporations

- ・Compliance with Industrial Safety and Health Act
- ・Experimental safety management and education

2. Stricter regulations in life science research

- ・Recombinant DNA experiments: Act on the Conservation and Sustainable Use of Biological Diversity through Regulations on the Use of Living Modified Organisms, etc.
- ・Experimental animals: Act on Welfare and Management of Animals, etc.
- ・Experiments handling pathogenic microorganisms: Act on the Prevention of Infectious Diseases and Medical Care for Patients with Infectious Diseases, etc.

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Request in the use of CRL facilities

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<Basic Rules>

1. Fill out the booking sheet of each instrument that you plan to use; Write the proposed time of use, your name and belonging, and the extension number.
2. Prepare any consumables, tools, chemicals, etc. for your own use.
3. Strictly follow the operating rules to use instruments and analyzers.
4. Fill out the user record sheet after every use (This is important for calculating the user charge).
5. Maintain the facility in a clean condition after use.

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Cautionary notice for booking equipment

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Follow the CRL operating rules to use equipment

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
Accident case 1

The rotor came off from the spindle and rolled around in the chamber. The rotor and spindle were severely damaged.

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Accident case 2

The aluminum screws of the tube caps were broken down by the centrifugal force at 50,000 rpm. Compatible tube caps must be used for ultracentrifugation.

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The bucket came off from the running rotor at 40,000 rpm and crashed onto the chamber wall at a velocity of over 1,600 km/h. Relative Centrifugal Field is 193,000 g.

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Tips for safety centrifuge operation

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- ・Central Research Laboratory Seminar (TBA)
- ・Technical Seminar (TBA)
- ・Intensive Course Seminar (held in September)
in "Basic Science Fundamentals & Multidisciplinary Seminars"

Participants in the above seminars can obtain credits for the PhD course lectures.

Details are indicated at the CRL website.
(<http://www.crl.shiga-med.ac.jp/home/seminar/seminar.html>)

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1. Please visit the CRL website for better usage of Research Equipment Section
(<http://www.crl.shiga-med.ac.jp/index.html>)
2. Good communication with CRL staffs is very helpful for both users and staffs.

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- ・ Users are expected to keep rules and regulations at CRL for obtaining accurate experimental data.
- ・ CRL staffs would greatly appreciate it if the users bring a copy of an article (PDF) to CRL when it is published.

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On site attendee

After this lecture, download the application form and submit it to CRL via e-mail

On line Zoom attendee

We will display a link for download site of two forms to be submitted ("survey form" and "application form for CRL User Registration") in the Zoom chat box during this course. By sending them using email, your registration will be completed.

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SDS (Safety Data Sheet)
化学物質や化学物質が含まれる原材料などを安全に取り扱うために必要な情報を記載したもの
SDS is a form containing data regarding the properties of a particular substance. It is intended to provide workers with procedures for handling or working with that substance in a safe manner.

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Introduction of Chemical Registration Information System

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Register information of all designated chemicals and high pressure gas cylinders in Chemical Registration Information System (CRIS FOREST)

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Update registered information in CRIS FOREST when designated chemicals and/or high pressure gas cylinders are used

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Introduction of Chemical Registration Information System (No.2)

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All work with formaldehyde must be performed in the fume hood of Morphological Research Joint Lab on the 4th floor. Formaldehyde Users should sign in a Formaldehyde Log after use.

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Transporting Liquid Nitrogen Container

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Emergency Procedures & Evacuation at CRL (in case of FIRE)

1. 助けを求める
Ask help without firefighting operations alone.
2. 防災監視室(内線2773)もしくはセンター職員に通報・連絡する(内線2300)
Contact immediately Security Control Center (Ext.2773) and/or the CRL office (Ext.2300).
3. 可能であれば初期消火、避難誘導、危険物の搬出に努める
Help extinguish the fire if it is at an early stage, provide evacuation instructions and bring out hazardous materials (e.g. alcohol).

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Evacuate from the fire zone if you think it dangerous.

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Emergency Procedures & Evacuation at CRL (in case of ACCIDENT)

1. まず自分の安全を確保し二次災害を防ぐ
Make sure you are safe and then prevent the secondary accident.
2. 怪我人がいる場合は適切な処置を行う
Help injured persons if it is possible.
3. センター職員に連絡する(内線2300)
Notify the condition to the CRL staffs (Ext.2300) immediately.
4. 後日、事故報告書を支援センターに提出
Submit an accident report to CRL.

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

How to use the recombinant DNA laboratory

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拡散防止措置(二種省令)
環境中への拡散を防止して行う使用等(大学等の実験室内)
Containment measures (Type 2 Use of LMOs)
Possible uses are those with containment measures to prevent to dispersal of LMOs in the environment (Uses in laboratories etc.)

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




P1レベル拡散防止措置(抜粋)
Containment Measures of P1 Level Laboratory


P1Aレベル拡散防止措置(抜粋)
Containment Measures of P1A(experiment with Animals) Level Laboratory

遺伝子組換え生物等が逃亡その他拡散しない構造の容器に入れること
When a LMO is taken out of the laboratory in process of an experiment, the LMO shall be put in a container of the structure that prevents it from leaking or other dispersion.

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P2レベル拡散防止措置(抜粋)
Containment Measures of P2 Level Laboratory

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
Application Process for Using Pathogenic Virus Lab

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Precautions of the recombinant DNA laboratory

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When using equipment for the first time, you must ask the CRL staffs how to operate it. Especially you must do so to use centrifuges.

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4) How to use experimental animals at CRL

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Need to complete following:
1. Completed the training program in animal experiment.
2. Pass the qualification exam of animal experiment.
3. Obtain the certification on animal experiment.

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Using animal laboratories in CRL Research Equipment Section

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Perform animal experiments under the Law Concerning the Protection and Control of Animals and the regulations of Research Center for Animal Life Science in SUMS.

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Comply the related law and regulations when performing animal experiments with X-ray and/or recombinant DNA.

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Approved animal laboratories on the 1st floor in CRL Research Equipment Section

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●エックス線照射装置等:
登録完了
Registration is completed for users
of X-ray irradiator

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Required items of education and training (4.5 hrs. ~)

- ① X-ray radiography work, ② Mechanism and safe handling of X-ray or gamma-ray irradiator,
- ③ Effects of ionizing radiation on the human body, ④ Laws and regulations

※Those who have already been registered as an X-ray radiography worker at the hospital may be exempted from items other than ① and ②.

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Required items of education and training (1 hr. ~)

- ① X-ray radiography work
- ② Mechanism and safe handling of X-ray or gamma-ray irradiator

学内登録申請(まとめ) Registration as a radiation worker in SUMS

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学内登録申請 Registration as a radiation worker in SUMS

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利用申込書の提出
Fill in a form for using X-ray generator

利用の承認
Acceptance for your use by CRL

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Equipment introduction #1 (Instruments on the 1st and 4th floors)

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- 大学内における実験・研究等に使用するヒト由来サンプルを一括して保管する施設
HSR is a facility where human-derived materials used for research activities within the University are securely and appropriately stored.
- 冷凍庫の温度は温度管理システムにより監視
HSR freezers are connected to temperature control systems.
- 超低温槽の異常時には利用者にメール
HSR registered users can get their freezers' temperature information from any area within the University. Also, they can be contacted by these alarm systems in the event the freezer is having a temperature critical issue.

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Laboratories and instruments on the 1st floors

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The device can be used for angiography, catheter positioning, and implantable transducer placement in small and medium-sized animals.

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Laboratories and instruments on the 1st floor

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
This equipment is used to irradiate cultured cells, microorganisms, and small animals such as mice and rats with large amounts of X-rays for the purpose of examining the effects of X-rays on living organisms.

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Laboratories and instruments on the 1st floors

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The system is a clinical wide-bore MRI available for MRI imaging of small, medium and large animals. Scanning experience includes crab-eating macaques, pigs, dogs, rabbits, rats, and mice (including transgenic mice).

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Laboratories and instruments on the 1st floor

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It is a device for estimating the structure of compounds by observing the phenomenon in which nuclei absorb specific electromagnetic waves (resonance phenomenon) corresponding to their respective chemical environments when nuclei are irradiated with electromagnetic waves from the outside in a given magnetic field.

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Laboratories and instruments on the 1st floors

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Wet Lab is a facility for animal experiments approved by Research Center for Animal Life Science and the Animal Experimentation Committee, and for surgical technique training of medical staff.

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Laboratories and instruments on the 1st floor

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This room is equipped with tools for producing labware and other equipment.

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Transmission Electron Microscope

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When the electron beam emerges from the specimen, the electron beam carries information about the structure of the specimen that is magnified by the objective lens system of the microscope.

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
The SEM produces images by probing the specimen with a focused electron beam that is scanned across a rectangular area of the specimen (raster scanning). Because the SEM images the surface of a sample rather than its interior, the electrons do not have to travel through the sample.

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Confocal laser scanning microscopy is a technique for obtaining high-resolution optical images with depth selectivity. The key feature of confocal microscopy is its ability to acquire in-focus images from selected depths, a process known as optical sectioning.

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
Dragonfly is high speed confocal microscope system. It combines deep imaging with a spinning disk system and high-speed scanning, and delivers outstanding multi-dimensional images, while significantly boosting productivity.

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
DeltaVision Elite is a fully integrated, flexible, widefield fluorescence microscope optimized for live cell imaging. It's efficient light path and optimized deconvolution algorithm combine to deliver superior image quality while maintaining cell viability.

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冷却CCDカメラ付倒立型蛍光・可視顕微鏡 Inverted Fluorescence Microscope with Cooled CCD Camera

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IX83 inverted microscope provides rapid and high-resolution imaging. It achieves time-lapse imaging of live cells with less phototoxicity and photobleaching.

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
Microscopy Sample Preparation Equipment

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This device automatically performs dehydration, de-alcoholization, and paraffin immersion of formalin-fixed tissues within a certain period of time.

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
This device is for embedding paraffin-immersed specimens in paraffin blocks.

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A microtome is a cutting tool used to produce extremely thin slices of material known as sections.

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A cryostat is a cutting tool used to produce extremely thin slices of frozen tissue known as sections.

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These are able to analyze many thousands of particles per second, in "real time" and, if configured as cell sorters, can actively separate and isolate particles with specified optical properties at similar rates. To analyze solid tissues, a single-cell suspension must first be prepared.

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Laboratories and instruments on the 4th floor

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Fluorescence activated cell sorting (FACS) utilizes flow cytometry to separate cells based on morphological parameters and the expression of multiple extracellular and intracellular proteins. This method allows multiparameter cell sorting and involves encapsulating cells into small liquid droplets which are selectively given electric charges and sorted by an external electric field.

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Laboratories and instruments on the 4th floor

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It is a cellular metabolism analyzer that enables non-invasive and highly sensitive serial measurement of the state of glycolysis and aerobic respiration by mitochondria, the main energy metabolism pathways in cells.

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Laboratories and instruments on the 4th floor

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It is used to culture cells under more physiological conditions by creating a stable environment with constant temperature, humidity, and CO₂ concentration, and is mainly used for the cultivation and growth of cells and microorganisms.

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Basic rules

- 予約が必要な機器は、予め予約表に記入する
A reservation is required for the use of the equipment at CRL.
Fill out the booking sheet of each instrument that you plan to use.
- 初めて利用する時は、必ず、利用方法、操作方法、注意事項を機器担当者に問い合わせる
When using equipment at CRL for the first time, any users should ask the CRL staffs how to operate it.
- 実験室に持ち込んだものは実験終了後持ち帰る
Do not leave your items in the working area, including waste materials (Keep the facility clean after use).

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8) Equipment introduction #2 (Instruments on the 2nd and 3th floor)

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Compatible to 6-384 well Plates for measuring absorption, fluorescence, and luminescence.

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Can measure DNA, RNA, and protein from 1–2 μ L of sample in seconds.

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
For protein detection on Western blot membranes.

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For quantitative analysis of gene expression.

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
For identification and quantification of up to 100 different molecules simultaneously in a single biomolecular assay.

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For tissue homogenization, cell lysis, DNA or chromatin shearing, and other sample preparations.

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Information such as concentration, size, and quality can be obtained as digital data by electrophoretic analysis of DNA, RNA, and protein.

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The device is capable of measuring the molecular weight of proteins, DNA, sugars, etc.

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Various rotors are available. For details, please visit our website.

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For detection of in vivo luminescence and fluorescence of small animals.

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Instruments on the 3rd floor

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Cold chambers for biochemical experiments and storage

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
For sample storage and biochemical experiments such as chromatography at 4°C.

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For disruption of chromatin, cells, bacteria, or tissue

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For homogenization of tissues uniformly in a short time.

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For drying samples without denature under freezing Temperature.

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For acquisition of anatomical, functional, physiological and molecular level data in vivo in real time.

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From April 2019, regular USB memory devices are not allowed to transfer the data in all computers in the CRL.
We recommend to use the USB memory device with encryption function.

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We recommend Kingston Technology's DataTraveler2000.

Equipment List (Research Equipment Section, Central Research Laboratory)

BIOCHEMISTRY	EQUIPMENT NAME	FEE (JPY)	LOCATION
Microplate Reader	TECAN/Infinite M PLEX	100/use	2F Biochemical Analysis Lab 1
Micro Mass Spectrometer	Thermo Fisher/NanoDrop One, NanoDrop ND-1000	—	2F Biochemical Analysis Lab 1
Recording Spectrophotometer	JASCO/V-570DS	—	2F MS Spectrometry Lab
Fourier Transform Infrared Spectrometer	SHIMADZU/IRPrestige-21	—	2F MS Spectrometry Lab
Fluorescence Spectrometer	HITACHI/F-2500	—	2F MS Spectrometry Lab
Circular Dichroism spectrometer	JASCO/J-1500	—	2F Biochemical Analysis Lab 1
Gel Imaging System	NIPPON Genetics/FAS5	—	2F Molecular Biology Analysis Lab
Nuclear Magnetic Resonance Spectroscopy	JEOL/JNM-ECZ400S, JNM-ECX400	10,00/h (Day Max 5,000)	1F NMR Lab
High Performance Liquid Chromatograph	SHIMADZU/HPLC system	200/use	2F Biochemical Analysis Lab 1
HPLC AKTA System	GE	200/use	2F MS Spectrometry Lab
Preparative Ultracentrifuge	Beckman/Optima L-90K	100/10min (5,000/More than 500min)	2F Centrifuge Room
Tabletop Ultracentrifuge	Beckman/Optima MAX-TL	100/10min (5,000/More than 500min)	2F Centrifuge Room
High Speed Refrigerated Centrifuge	Beckman/HP-25·HP-26	200/use	2F Centrifuge Room
Hybrid High Speed Refrigerated Centrifuge	KUBOTA/6200	—	2F Centrifuge Room
Low Speed Centrifuge	Himac/CF7D2	—	2F Molecular Biology Analysis Lab
High Speed Refrigerated Centrifuge (P2 Lab)	Beckman/Avanti J-E	200/use	2F Gene Engineering Lab P2
High Speed Refrigerated Microcentrifuge	TOMY/MX-305	—	2F Gene Engineering Lab P2
	HITACHI/CF16RN	—	3F Centrifuge Room
Plate Centrifuge	KUBOTA/PlateSpin II	—	2F Molecular Biology Analysis Lab
Multipurpose Centrifuge	Beckman/GS-15R	—	2F Gene Engineering Lab P2
Low Speed Refrigerated Centrifuge	Himac/CF7D2	—	2F MS Spectrometry Lab
Centrifugal Concentrator/Speedvac Concentrator	Thermo Fisher/SPD1010	200/use	2F Centrifuge Room
Deep Freezer	NIHON FREEZER/CLN-52U	6,000/y	2F Gene Engineering Lab P2
Shaking Incubator	Sanki Seiki/SCS-R·SCS-12R, TAITEC/BR-30L,BR-30LF,BR-43FL	—	2F Gene Engineering Lab P2
Vacuum Freeze Drying Equipment	LABCONCO/LL-1	200/d	3F Cell homogenization Room
Ultrasonic Homogenizer	BRANSON/Sonifier Model 450-Advanced	—	3F Cell homogenization Room
Cap Horn Homogenizer	MISONIX/ASTRASON MODEL XL2020	—	3F Cell homogenization Room
Polytron Homogenizer	KINEMATICA/MODEL PT1200E	—	3F Cell homogenization Room
Chromatochamber	NIHON FREEZER/MC-30EF3	500/rack/m	2F Frozen Storage Room
Electronic Balance	SHIMADZU/AUW320	—	2F Molecular Biology Analysis Lab

BIOCHEMISTRY	EQUIPMENT NAME	FEE (JPY)	LOCATION
Gas Chromatograph (EI)	SHIMADZU/GCMS-QP5050A	100/use	2F MS Spectrometry Lab
Time-of-Flight Mass Spectrometer	SHIMADZU/MALDI-TOF/MS AXIMA Confidence	1,000/use	2F MS Spectrometry Lab
Inductively Coupled Plasma Mass Spectrometer	Agilent/7700e ICP-MS	Basic charge 800 +1,400/h	2F MS Spectrometry Lab
X-ray Diffractometer	Philips/MPD-1880	1,000/use	1F X-ray Irradiation Lab
DNA Sequencer	ABI/PRISM 3130XL GeneticAnalyzer HITACHI/DS3000	200/sample	3F Gene Analysis Lab
Acoustic Solubilizer	Covaris/S220	100/use	2F Molecular Biology Analysis Lab
Real Time Quantitative PCR (LightCycler)	Roche/LightCycler480 System II	100/use	2F Molecular Biology Analysis Lab
Real Time Quantitative PCR (Plate)		7,600/set	2F Molecular Biology Analysis Lab
Thermal Cycler	ABI/GeneAmp PCR System 2720, Roche/LightCycler	—	2F Molecular Biology Analysis Lab
	Thermo Fisher/MiniAmp Plus	—	3F Molecular Biology Analysis Lab
	Applied Biosystems/GeneAmp PCR System 9600	—	3F Analysis Room
Lumino Image Analyzer	FUJIFILM/LAS-4000	100/use	2F Molecular Biology Analysis Lab
Chemiluminescence Imaging System	Vilber-Lourmat/FUSION	100/use	2F Molecular Biology Analysis Lab
Electrophoresis Apparatus	Pharmacia/IPG Phor	500/use	2F MS Spectrometry Lab
Bioanalyzer	Agilent Technologies/Agilent 2100 Bioanalyzer Limited	100/use	2F Molecular Biology Analysis Lab
Suspension Array System	Bio-Rad/Bio-Plex200	3,000/use	2F Biochemical Analysis Lab 1
MOE (Molecular Operating Environment)	Chemical Computing Group/MOE	—	3F Analysis Room
Heating Stirring Module	ThermoFisher/Reacti-Therm	200/h	3F Analysis Room
Sample Concentrator		200/h	3F Analysis Room

MORPHOLOGY	EQUIPMENT NAME	FEE (JPY)	LOCATION
Transmission Electron Microscope	HITACHI/H-7500	500/h + film cost (Market Price)	4F Electron Microscope Lab
Field Emission Scanning Electron Microscope	JEOL/JSM-7505FA	2,000/h	4F Electron Microscope Lab
Scanning Electron Microscope	JEOL/JSM-6010LA	2,000/h	4F Electron Microscope Lab
Ultramicrotome	Reichert/ULTRACUTE	500/use	4F Microscopy Sample Preparation Room 1
Cooled CCD Camera	NIPPON ROPER/CoolSNAP-HQ	300/use	4F Electron Microscope Lab
Color CCD Camera	NIPPON ROPER/Retiga2000R Nikon/ECLIPSE Ni-E	100/use	4F Cell Culture Room 4F Image Processing Room
High-resolution Fluorescence Microscope	GE/DeltaVision Elite	200/use	4F Microscope Lab 1
All-in-One Fluorescence Microscope	KEYENCE/BIOREVO BZ-9000	—	4F Cell Culture Lab
Confocal Laser Scanning Microscope	OLYMPUS/FV1000-D	1,000/use	4F Microscope Lab 2
	LEICA/TCS SP8 X	1,000/use	4F Confocal Laser Scanning Microscope Lab 3
Inverted Fluorescence Microscope	OLYMPUS/IX83	100/use	4F Microscopy Sample Preparation Room 2

MORPHOLOGY	EQUIPMENT NAME	FEE (JPY)	LOCATION
Live Cell Imaging Microscope	Andor/Dragonfly200	300/h (Day Max 5,000)	4F Confocal Laser Scanning Microscope Lab 1
Microtome	LEICA/SM2010R	100/use	4F Microscopy Sample Preparation Room 2
Tissue Processor	SAKURA/Tissue-Tek VIP 6-J0	1,000/use	4F Microscopy Sample Preparation Room 2
Tissue Embedding Console	SAKURA/Tissue-Tek 6	100/piece	4F Microscopy Sample Preparation Room 2
Cryostat	LEICA/CM3050 S	500/use	4F Microscopy Sample Preparation Room 2
3D/4D Image Analysis and Visualization Software	BITPLANE/Imaris	100/use	4F Image Processing Room
Sample Preparation	Paraffin embedding	500/use	4F Microscopy Sample Preparation Room 1
	Sectioning	300/slide	4F Microscopy Sample Preparation Room 1
	H&E Staining	300/slide	4F Microscopy Sample Preparation Room 1

CELL CULTIVATION	EQUIPMENT NAME	FEE (JPY)	LOCATION
Electroporator	LONZA/4D-Nucleofector, amaxa/Nucleofector	100/use	4F Cell Culture Lab
CO2 Incubator	YAMATO/IP400, SANYO/MCD-175, ESPEC/BNA-111, HIRASAWA/CPD-2701	45,000/y + CO2 cost	4F Cell Culture Lab
CO2		6,000/cylinder	4F Cell Culture Lab
Clean Bench (Laminar Flow Cabinet)	ESCO/Airstream PRO ESC-BCB-4A7, Showa Kagaku/S-1801WBV, etc.	—	4F Cell Culture Lab
Flexercell	FLEXERCELL/FLEX-I	2,000/use	4F Cell Culture Lab
Microinjection System	Eppendorf/InjectMan NI2, OLYMPUS/IX-73	100/use	4F Cell Engineering Lab 3
Cell Analyzer	BD/FACSCalibur・FACSCanto II	1,500/use	4F Cell Engineering Lab 2
	BD/LSRFortessa	1,500/use	4F Cell Engineering Lab 4
Cell Sorter	BD/FACSAria	5,000/use (Analysis)1,500/use	4F Cell Engineering Lab 1
	BD/FACSAria Fusion	5,000/use (Analysis)1,500/use	4F Cell Engineering Lab 4
Magnetic Cell Sorter	Miltenyi Biotec/autoMACS Pro Separator	3,500/use (No consumables)1,100/use	4F Cell Engineering Lab 1
Cell Dynamic Analysis System	GE/EZ-TAXIScan	unfixed	4F Cell Culture Lab
Extracellular Flux Analyzer	Seahorse/XFe24	100/use	4F Cell Culture Lab
Deep Freezer (Frozen Strage Room) Container Basket	REVCO/ULT-1490-3J-D30, etc.	(Container)6,000/y (Drawer Case)4,000/y	3F Frozen Strage Room
Large Liquid Nitrogen Sample Storage Container	CHART/MVE Heco1536P-190AF・GB	550/m	4F Workroom
Biosafety Cabinet	ESCO/ClassII Type A2	—	4F Cell Engineering Lab 4
	Thermo/1387	—	4F Pathogenic Virus Lab
Preparative Ultracentrifuge	Beckman/Optima L-90K	100/10 min (5,000/More than 500min)	4F Pathogenic Virus Lab
High Speed Refrigerated Centrifuge	Beckman/Avanti J-E	200/use	4F Pathogenic Virus Lab
Tabletop Centrifuge	KUBOTA/KS-5000	—	4F Pathogenic Virus Lab
Tabletop Refrigerated Centrifuge	Beckman/Allegrax-30R	—	4F Cell Culture Lab


CELL CULTIVATION	EQUIPMENT NAME	FEE (JPY)	LOCATION
Ethylene Oxide Gas Sterilizer	ELK CORPORATION/ SA-360ECO	3,000/use	4F Frozen Strage Room
X-ray Irradiator	AcroBio/RS320	1,000/use	Research Center for Animal Life Science

PHYSIOLOGY / ANIMAL	EQUIPMENT NAME	FEE (JPY)	LOCATION
Microdialysis	EICOM/HTEC-500	1,000/use	2F Biochemical Analysis Lab 1
Animal Blood Gas Analyzer	ABAXIS/VetScan VS2	—	3F Biological Sample Analysis Lab
Automatic Blood Cell Counter	NIHON KOHDEN/ Celltac α MEK-6550	200/sample	3F Biological Sample Analysis Lab
Oxidation Redox Analyzer	WISMERLL/REDOXLIBRA	100/use	3F Biological Sample Analysis Lab
Laboratory Animal X-ray Fluoroscope	TOSHIBA Medical/Plessart50	1,000/use	1F Animal X-ray Fluoroscopy Lab
Langendorff System	ADInstruments	1,000/use	4F Cell Engineering Lab 3
High-resolution Ultrasonic Echo Imaging System for Small Animals	VisualSonics/Vevo2100	250/15min (Day Max 5,000)	3F Small Animal Imaging Analysis Lab
In vivo Imaging System for Small Animals	Vilber-Lourmat/NEWTON 7.0	200/h	3F Small Animal Imaging Analysis Lab
3-Tesla MRI Scanner	SIEMENS Healthineers/ MAGNETOM Verio dot 3T	400,000/y 5,000/d, 1,500/h	1F 3T MRI Lab
4.7-Tesla MRI Scanner	BRUKER/BioSpec47/40USR	400,000/y 5,000/d, 1,500/h	Research Center for Animal Life Science


ADDITIONAL SERVICE	EQUIPMENT NAME	FEE (JPY)	LOCATION
Elix Water	MILLIPORE/Milli-Q Integral 5	10/ L	2F Biochemical Analysis Lab 4F Workroom
Milli-Q Water	MILLIPORE/Milli-Q Integral 5	80/ L	2F Biochemical Analysis Lab 4F Workroom
Ice Machine	Scotsman/F0522B Hoshizaki/FM-340AK-SA	—	2F Biochemical Analysis Lab 4F Workroom
Liquid Nitrogen	100 liter, 50 liter tank	400/L (Market Price)	4F Workroom
Oxygen Gas	purity 99.5%, 4.0Kg/cm ²	3,630/cylinder (Market Price)	1F Bombe Room 1F Wet Lab
Compressed Air		5,500/cylinder (Market Price)	1F Wet Lab
Human Sample Storing	REVCO, NIHON FREEZER	2,000/freezer/m	1F Human Sample Resource Room
Space Charge		2,000/3.3m ² /m	
Equipment Rental		200/w	

WET LAB	EQUIPMENT NAME	FEE (JPY)	LOCATION
Electric Knife	ellman/SURGITRON Dual EMC 90	—	1F Wet Lab
Wet Lab Space Rental		2,000/3.3m ² /m	1F Wet Lab

1F·2F of Central Research Laboratory

3F·4F of Central Research Laboratory

● Access Management System

Contact Number & E-mail Address

Research Equipment Section of Central Research Laboratory

○ Director of CRL

Job title	Name	Ext.	E-mail
Professor	Itoh, Yasushi	2171	yasushii@belle.shiga-med.ac.jp

○ Faculty of CRL

Job title	Name	Ext.	E-mail
Associate Professor	ASAHINA, Kinji	2301	asahina@belle.shiga-med.ac.jp

○ Technical Support Staff

Job title	Name	Ext.	E-mail
Technical Staff	MORI, Yasuhiro	2307	moriyasu@belle.shiga-med.ac.jp
Technical Staff	TERADO, Tokio	2306	terado@belle.shiga-med.ac.jp
Technical Staff	YAMAMOTO, Takefumi	2304	tyama@belle.shiga-med.ac.jp
Technical Staff	URUSHIYAMA, Noboru	2302	urushi@belle.shiga-med.ac.jp

○ Administrative Staff

Job title	Name	Ext.	E-mail
Administrative Staff	YAMAKAWA, Nobuko	2300	nyama@belle.shiga-med.ac.jp

○ Central Research Laboratory

Location	Ext.	E-mail
2F Central Research Laboratory office	2300	hqcr1@belle.shiga-med.ac.jp

More information on CRL is available at the CRL web site (<http://wwwcrl.shiga-med.ac.jp/index.html>) where you can get a list of installed equipment, online edition of the operation manuals, etc. and download various request/application forms.


Application Procedures for New User Registration

Those who wish to use the facilities and equipment of Central Research Laboratory (CRL) must be registered users. To be the user, they are required to take the Guidance for New Users managed by CRL. Depending on the nature of the experiments they plan to conduct, other application procedures may be required in addition to the user registration.

They must complete the necessary procedures according to the section they wish to use and the type of experiments they are planning to conduct.

Once you have completed the registration ···

Each piece of equipment at the CRL has its own rules on how to reserve and use it. If you use equipment for the first time, please be sure to ask the staff member in charge of the equipment how to use it before using it, rather than having your lab members teach you how to use it.

Items that Require Application to the University

- ✓ Recombinant DNA Experiments: Application to the SUMS Genetic Recombination Experimentation Safety Committee (Research Promotion Division)
- ✓ Experiments Handling pathogenic microorganisms: Application to the SUMS Biosafety Committee (Research Promotion Division)
- ✓ Experimental Animals: Application to the SUMS Animal Experimentation Committee (Research Center for Animal Life Science)
- ✓ Experiments with Animal X-ray Fluoroscope: Registration as an X-ray photographing worker (Research Promotion Division)
- ✓ Experiments with X-ray Irradiator: Registration as an X-ray worker (Research Promotion Division)
- ✓ Human Subjects Research: Application to the SUMS Ethics Committee

Please follow the instructions below when completing “Application Form for CRL User Registration”

Those who will conduct the following experiments at CRL must enter the approval number and other necessary information on the Form and attach the necessary documents, if necessary:

- (1) **Recombinant DNA experiments**: The approval number must be entered. If you plan to use the Gene Engineering Lab, please contact the staff (Ext.2306) to register your ID card information to the entry/exit managing system.
- (2) **Experiments handling pathogenic microorganisms**: The approval number must be entered and the prescribed documents must be submitted. If you plan to use the Pathogenic Virus Lab and/or the Cell Engineering Lab #4, please contact the staff (Ext.2307) to register your ID card information to the entry/exit managing system.
- (3) **Animal experiments**: The approval number must be entered.
- (4) **Experiments using X-ray generators (Animal X-ray Fluoroscope, X-ray Irradiator, X-ray Diffractometer)**: X-ray generator users must take the X-ray training course.
- (5) **Human Sample Resource Room**: The prescribed form must be submitted. Persons using the HSR must register their ID card information to the entry/exit managing system prior to using the room.

【Caution】 Matters to be entered, documents to be submitted

Items	Approval number	Entry/exit managing system	Document to submit	Documents to submit	Staff (Ext.)
Recombinant DNA experiments	●	●			Terado (2306)
Handling pathogenic microorganisms	●	●	●	①CRL Application Form for Handling Pathogens, ②Application for Handling Pathogens (Form 1 or 2)	Mori (2307)
Experimental Animals	●				Terado (2306)
Animal X-ray Fluoroscope		●			Yamamoto (2304)
Wet Lab	●				Terado (2306)
Human Sample Resource Room		●	●	HSR Access Application	Yamamoto (2304)

EXAMPLE Take extra care not to make mistake or omit anything

Central Research Laboratory (CRL), Research Equipment Section
2023 New User Registration Form

Date: 2023/00/00

I will abide by rules and regulations established by CRL.

Affiliation	Department/Center of 00
Position	Graduate student
Name	SHIGA Taro
Staff (Student) ID number	88888888
Phone number	2300
e-mail address	hqcr1 @belle.shiga-med.ac.jp
Supervisor's name	

Brief description of your research projects:

☒ of equipment that you plan to use

- | | | | |
|---|--|--|---|
| <input type="checkbox"/> Spectrophotometer | <input type="checkbox"/> Bio-Plex suspension array | <input type="checkbox"/> Live imaging microscope | <input type="checkbox"/> Heart perfusion system |
| <input type="checkbox"/> Mass spectrometer | <input type="checkbox"/> MOE program | <input type="checkbox"/> Fluorescence microscope | |
| <input type="checkbox"/> Nuclear magnetic resonance | <input type="checkbox"/> Chemilumi imaging | <input type="checkbox"/> CO ₂ incubator | Take the X-ray training course |
| <input type="checkbox"/> Temperature control centrifuge | <input type="checkbox"/> PCR/qPCR | <input type="checkbox"/> MACS | <input type="checkbox"/> VetScan |
| <input type="checkbox"/> Ultracentrifuge | <input type="checkbox"/> Gel imaging | <input type="checkbox"/> Flow cytometer | <input type="checkbox"/> X ray fluoroscope |
| <input type="checkbox"/> Lyophilizer | <input type="checkbox"/> Chromatography chamber | <input type="checkbox"/> FACS | <input type="checkbox"/> X ray generator |
| <input type="checkbox"/> Vaccume concentrator | <input type="checkbox"/> Electronmicroscope | <input type="checkbox"/> EZ-TAXIScan | <input type="checkbox"/> MRI |
| <input type="checkbox"/> Micro-dialysis analyzer | <input type="checkbox"/> EM sample preparation | <input type="checkbox"/> Nucleofector | <input type="checkbox"/> Millipore water |
| <input type="checkbox"/> Homogenizer | <input type="checkbox"/> Light microscope | <input type="checkbox"/> Flux analyzer | <input type="checkbox"/> Liquid nitrogen |
| <input type="checkbox"/> DNA sequencer | <input type="checkbox"/> Microscope sample prep | <input type="checkbox"/> Deep freezer | <input type="checkbox"/> Oxygen gas |
| <input type="checkbox"/> Bioanalyzer | <input type="checkbox"/> Image analyzer | <input type="checkbox"/> Liquid nitrogen dewar | <input type="checkbox"/> Handicraft tools |
| <input type="checkbox"/> Covaris solubilizer | <input type="checkbox"/> Laser microscope | <input type="checkbox"/> Autoclave | <input type="checkbox"/> Human sample storage |

Please check boxes if your experiment at CRL is involved in recombinant DNA, pathogenic viruses and microorganisms, or animal use.

※ New users will need to submit an additional application for entering these room by a card key.

<input checked="" type="checkbox"/> Recombinant DNA (<input type="checkbox"/> New • <input type="checkbox"/> Renewal)	Approved number : 04-1		Supervisor's name : (Affiliation :)	
	Approved room	<input type="checkbox"/> Genetic Engineering Lab (※) <input type="checkbox"/> X-ray Irradiation Lab <input type="checkbox"/> 3T MRI Lab <input type="checkbox"/> Cell Culture Lab <input type="checkbox"/> Animal Imaging Lab <input type="checkbox"/> Laser Microscope Lab 1 <input type="checkbox"/> Cell Engineering Lab 1 <input type="checkbox"/> Cell Engineering Lab 2 <input type="checkbox"/> Cell Engineering Lab 3 <input type="checkbox"/> Cell Engineering Lab 4 (※)		
<input checked="" type="checkbox"/> Pathogenic viruses Contact the staff (Ext.2307) and follow the prescribed procedures	Approved number : 04-1		Supervisor's name : (Affiliation :)	
	Approved room	<input type="checkbox"/> Pathogenic Virus Lab (※) <input type="checkbox"/> Cell Engineering Lab 4 (※)		
<input checked="" type="checkbox"/> Animal experiments (<input type="checkbox"/> New • <input type="checkbox"/> Renewal)	Approved number : 2022-1		Supervisor's name : (Affiliation :)	
	Approved room	<input type="checkbox"/> X-ray Fluoroscopy Lab (※) <input type="checkbox"/> X-ray Irradiation Lab <input type="checkbox"/> Wet Lab <input type="checkbox"/> 3T MRI Lab <input type="checkbox"/> Animal Imaging Lab <input type="checkbox"/> Laser Microscope Lab 1 <input type="checkbox"/> Cell Engineering Lab 3		

[NOTE] Any person who is engaged in the operation of recombinant DNA experiments, experiments handling pathogenic microorganisms, animal experiments must comply with the regulations in life science. research.