



Fusion(70u) 488, 561, 640, 405,355 (3-4-3-6-3)(BD).

Cytometer:	FACSAriaIII	User:	Administrator
Cytometer Name:	FACSAriaIII	Institution:	N/A
Serial Number:	1	Software:	BD FACSDiva 9.0.1
Input Device:	Manual	Date:	9/10/2021 5:01:23 PM
Sheath Pressure:	70.00		
Nozzle Size:	70		
Window Extension:	2.00		


Blue Laser
(488nm)
FSC




561 Yellow-Green Laser
(561nm)




Red Laser
(640nm)



Violet Laser
(405nm)



355 UV Laser
(355nm)



Comments: